

seca 354

Digital baby scale also
converts to a flat scale
for children



- Baby scale and floor scale for children in one.
- Weighing tray easy to detach from base.
- Breast-Milk-Intake-Function.
- Robust yet lightweight scale is simple to transport.

seca 354:

Digital baby scale also converts to a flat scale for children.

Wherever life begins in the world, whether it be at home or in surgery: reliable weight control is decisive for a healthy start in life. The seca 354 is an all-around scale that has been developed on the basis of practical experience. It can be used not only as a baby scale but also as a floor scale for small children.



Baby scale and floor scale for children in one.

Two scales in one! The seca 354 baby scale has a comfortable weighing tray and a fine graduation of 0.5 oz. In no time at all, it can be converted into a floor scale for children with a capacity of 44 lbs. The weighing tray and the base which serves as the floor scale are securely locked together, thus ensuring that the baby can be weighed safely and securely.

Very safe due to secure attachment of tray.

Pressing a button unfastens the mechanism attaching the tray to the base. When locked, this secure attachment ensures that the baby scale is easy to lift and safe to transport.

Easy to operate due to self-explanatory button functions.

The seca 354 is easy to operate due to its intuitive operating concept with self-explanatory buttons. The auto-HOLD function enables the baby to be taken care of before one has to read off the result. With the TARE function, the weight of a blanket or soft mat can be reliably deducted. A further advantage of the scale is its Breast-Milk-Intake-Function which determines how much milk the baby has consumed. The robust but lightweight scale is easy to transport and the user friendly buttons make it easy to operate.



Mobile and robust for transport.

Its compact, robust construction predestines the seca 354 as the scale for all those who are on-the-go and attending to the health of their youngest patients. The power-saving battery operation measures up to 20,000 weighings with just one set of batteries, making the scale independent of electrical power supply.



seca 413



seca 414



The weighing tray can be removed quickly and easily at the touch of a button.



seca 354

Technical Data

- Capacity: 44 lbs / 20 kg
- Graduation: 0.2 oz < 22 lbs > 0.5 oz / 5 g < 10 kg > 10 g
- Dimensions (WxHxD): 21.7 x 6.1 x 13.1" / 552 x 156 x 332 mm
- Dimensions/tray (WxHxD): 20.7 x 3.1 x 9.8" / 525 x 80 x 250 mm
- Dimensions/platform (WxHxD): 11 x 2.6 x 9.3" / 280 x 65 x 235 mm
- Weight: 5 lbs / 2.3 kg
- Power supply: Batteries
- Functions: TARE, BMIF, auto-HOLD, lbs/kg switch-over, automatic switch-off
- Optional: measuring board seca 417, carrying case seca 414 (for combination with measuring board) and seca 413