

seca 684

Digital multifunctional scale



- Safe weighing of standing, sitting or wheelchair patients.
- Large, easy-to-access platform.
- Wheel barrier provides extra safety.
- Drive-up ramp seca 470 included.
- TARE, auto-HOLD and damping function.
- Integrated transport castors allow for easy mobility.

seca 684:

Weighing without restrictions.

The multifunction scale seca 684 is seca's jack-of-all-trades. Thanks to its many functions, this scale is ideal for use in hospitals, general medicine, nephrology, youth medicine as well as in nursing homes. We have addressed the many needs of daily medical life through sophisticated and functional design. The result? A scale that makes your work easier, everyday.



Safe weighing of standing, sitting or wheelchair patients.

On this scale patients can stand aided by a stable handrail, sit on the integrated folding seat, or remain in a wheelchair on the large platform (31.5 x 33.5"). When standing, a stable handrail provides a safe hold for unstable patients. When sitting, a comfortable seat for minimally ambulatory patients can be raised and lowered with one hand. Special attention was given to making the scale easy to access, and a stop barrier for the front wheels of the wheelchair provides additional safety.

TARE and auto-HOLD function.

TARE and pre-TARE functions enable reliable net weight determination in spite of added weight, for example by deducting the known weight of a wheelchair. The practical auto-HOLD function shows the weight even after the patient has left the scale. The weight is always easy for patient and caregiver to read from the swiveling display mounted at convenient hip level.

A fold-out seat with integrated damping system?

There is such a thing – from seca of course. The multifunctional scale seca 684 is equipped with a stop barrier for wheelchairs and another detail that adds to overall safety. A spring in the mechanism for the integrated seat keeps the seat securely in place until it is released for use.

seca 684

Technical data

- Capacity: 800 lbs / 360 kg
- Graduation: 0.1 lbs / 50 g
- Dimensions (WxHxD): 35.8 x 43.3 x 41.9" / 910 x 1,100 x 1,064 mm
- Dimensions platform (WxHxD): 29.5 x 2 x 32.3" / 760 x 50 x 820 mm
- Weight: 130 lbs / 59.5 kg
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, kg/lbs switch-over, damping, SEND/PRINT
- Power supply: Power adapter
- Data transmission: seca 360° wireless technology
- Optional: RS232 adapter set seca 460,
- Including ramp seca 470
- System compatible with: 360° wireless printer advanced seca 466, 360° wireless printer seca 465, software seca analytics 105 and seca emr flash 101, 360° wireless USB adapter seca 456, stadiometers seca 274 and seca 264

Useful details for perfect handling.

The non-slip rubber coating and handrail add safety, while integrated transport castors make the scale easy to move.



You can stand or sit when you are being weighed.



Thanks to the integrated ramp, a wheelchair rolls easily and comfortably onto the platform.

