



WE CARE LIKE FAMILY™

# TRANSPORT SHARPS CONTAINERS



## Sharps Shaft

Reorder No. 4630

- Designed to help safely dispose of used syringes and sharps
- Snap-on, leak-resistant cap provides temporary closure between uses
- Puncture-resistant plastic helps prevent accidental needlesticks
- Secure locking mechanism allows for safe final disposal
- Lightweight and compact; ideal for emergency kits
- Translucent container provides view of contents to avoid overfilling
- Measures 6.37" x 1.47" x 3.46"
- 24/cs

## Sharps Crate

Reorder No. 4631

- Designed to safely dispose of used medical needles and other sharps
- Disposable crate is crafted from puncture-resistant plastic for improved safety
- Compact and lightweight for easy portability in travel bags and emergency kits
- Suitable for use by law enforcement agencies for the safe containment of hazardous evidence
- Measures 2.03" x 7.2" x 1.9"; Volume: 0.5 qt. (0.473 liters)
- 40/cs

