smiths medical medex

Pressure Monitoring

Product Catalog



Pressure Infusor Bags



Pressure Transducers



Closed Blood Sampling



Accessories and Components



Pressure Monitoring Cables



The Right Choice.

smiths medical bringing technology to life

	PRESSURE INFUSOR BAGS	4
	Clear-Cuff®	4
	C-Fusor® Pressure Infusors	4
	PRESSURE TRANSDUCER	7
	DISPOSABLE TRANSDUCER SYSTEMS	8
	TranStar® Pressure Transducer Single Kits	8
	TranStar® Pressure Transducer Double Kits	9
	TranStar® Pressure Transducer Triple Kits	9
	TranStar® Pressure Monitoring Kits	10
	REUSABLE TRANSDUCER SYSTEMS	11
	LogiCal® Pressure Transducer Single Kits	11
	LogiCal® Pressure Transducer Double Kits	12
	LogiCal® Pressure Transducer Triple Kits	12
	LogiCal® Pressure Transducer Cartridges	13
	LogiCal® Pressure Transducer System Component	13
	CLOSED BLOOD SAMPLING SYSTEMS	15
	Add-On Closed Blood Sampling	16
	TranStar® Kids Kit™ Neonatal Closed Blood Sampling System	18
	TranStar® HemoDraw® Plus Single Line Closed Blood Sampling System	19
	TranStar® HemoDraw® Plus Double Line Closed Blood Sampling System	
	TranStar® HemoDraw® Plus Triple Line Closed Blood Sampling System	
	LogiCal® Kids Kit™ Neonatal Closed Blood Sampling System	
	LogiCal® HemoDraw® Plus Single Line Closed Blood Sampling System	
	LogiCal® HemoDraw® Plus Double Line Closed Blood Sampling System	
	LogiCal® HemoDraw® Plus Triple Line Closed Blood Sampling System	
	ACCESSORIES AND COMPONENTS	
	TranStar® Mounting Accessories	
	LogiCal® Mounting Accessories	
	Blood Sampling Access Devices	
	HemoDraw® Plus Mounting Accessories	
	Flush Devices	
	Pressure Monitoring Tubing with Male/Female Luer Lock	
	Pressure Monitoring Tubing with Male/Male Luer Lock	
	Pressure Monitoring Tubing with Male/Female Luer Lock (Red)	
	Pressure Monitoring Tubing with Male/Female Luer Lock (Low Volume)	
	Non-Vented Administration Sets	
	Pressure Monitoring Components	
ì	Manometers	
	PRESSURE MONITORING CABLES	
	Pressure Monitoring Cables	
	INDEX	37






PRESSURE INFUSORS

Clear-Cuff®
C-Fusor® Pressure Infusors

Clear-Cuff®

ORDERING INFORMATION


Clear-Cuff®		
Code	Description	Units per Case
MX4705	Clear-Cuff® 500 mL Reusable Pressure	10
MX4710	Clear-Cuff® 1000 mL Reusable Pressure	10



C-Fusor® Pressure Infusors

ORDERING INFORMATION

C-Fusor® Disposable Dome			
Code	Description	Units per Case	
MX4805	C-Fusor® 500 mL Pressure Infusor	1	
MX4810	C-Fusor® 1000 mL Pressure Infusor	1	
MX4830	C-Fusor® Extension 3000 mL	1	
MX1821B	C-Fusor® Replacement Squeeze Bulb Assembly	1	
MX1821G	C-Fusor® Replacement Pressure Gauge	1	



p. 4 Pressure Monitoring



PRESSURE TRANSDUCERS

Disposable Transducer Systems



TranStar® Pressure Transducer Single Kits......8



LogiCal® Pressure Transducer Single Kits	1
LogiCal® Pressure Transducer Double Kits	12
LogiCal® Pressure Transducer Triple Kits	12
LogiCal® Pressure Transducer Cartridges	13
LogiCal® Pressure Transducer System Component	13

TranStar® Pressure Transducer Single Kits

ORDERING INFORMATION


Code	Description	Units per Case
4X9501T	TranStar® 12" (29 cm) Patient Mount Monitoring Kit	10
	12"	
MX9504T	TranStar® 60" (152 cm) Single Monitoring Kit	10
	48" 12"	
MX9505T	TranStar® 84" (213 cm) Single Monitoring Kit	10
	72"	
1X9510T	TranStar® 72" (183 cm) Single Line (without Administration Set), may be used with MX9513T	10
	60" 12"	
MX9530T	Neonatal TranStar® Monitoring Kit, 12", 30 cc/hr, Flush Device with Orange Insert (Neonatal)	10

p. 8 Pressure Monitoring

TranStar® Pressure Transducer Double Kits

ORDERING INFORMATION


TranStar® Pressure Transducer Double Kits		
Code	Description	Units per Case
MX9502T	TranStar ® 72" (183 cm) Double Monitoring Kit	10




TranStar® Pressure Transducer Triple Kits

ORDERING INFORMATION

TranStar®	Pressure Transducer Triple Kits	
Code	Description	Units per Case
MX9506T	TranStar® 72" (183 cm) Triple Monitoring Kit	10



MX9513T TranStar® 72" (183 cm) Trifurcated Add-On Kit, may be used with MX9510T 10



TranStar® Pressure Transducers

ORDERING INFORMATION


Code	Description	Units per Case
MX950	12" Cable	10
	A SESSION OF THE PROPERTY OF T	
4X950L	60" Cable	5
	60" 832	
MX950XSC	One Stopcock	10
	D A B A B A B A B A B A B A B A B A B A	
1X950X2SC	Two Stopcocks	10
	PRESSS	
1X950SCFT	With one Stopcock and 3 cc/hr Flush	10
	PRESS PRESS	
ИX950SCNT	With one Stopcock and 30 cc/hr Flush with Orange Insert (Neonatal)	10
	Press S	

p. 10 Pressure Monitoring

LogiCal® Pressure Transducer Single Kits

ORDERING INFORMATION

LogiCal® P	ressure Transducer Single Kits	
Code	Description	Units per Case
MX9601A	LogiCal® 12" (29 cm) Patient Mount Monitoring Kit	10



MX9604A LogiCal® 60" [152 cm] Single Monitoring Kit



MX9605A LogiCal® 84" (213 cm) Single Monitoring Kit



MX9630A Neonatal LogiCal® 12" Monitoring Kit, 30 cc/hr, Flush Device with Orange Insert (Neonatal) 10



LogiCal® Pressure Transducer Double Kits

ORDERING INFORMATION


LogiCal® Pressure Transducer Double Kits		
Code	Description	Units per Case
MX9602A	LogiCal® 84" (213 cm) Double Monitoring Kit	10



LogiCal® Pressure Transducer Triple Kits

ORDERING INFORMATION

LogiCal® Pressure Transducer Triple Kit		
Code	Description	Units per Case
MX9606A	LogiCal® 72" (183 cm) Triple Monitoring Kit	10



p. 12 Pressure Monitoring

LogiCal® Pressure Transducer Cartridges

ORDERING INFORMATION

50
25
50
10

LogiCal® Pressure Transducer System Component

ORDERING INFORMATION

LogiCal ® Reusable Transducer		
Code	Description	Units per Case
MX960	LogiCal® Reusable Transducer Mounting Plate	5







CLOSED BLOOD SAMPLING SYSTEMS


Add-On Closed Blood Sampling	. 16
TranStar® Kids Kit™ Neonatal Closed Blood Sampling System	. 18
TranStar® HemoDraw® Plus Single Line Closed Blood Sampling System	. 19
TranStar® HemoDraw® Plus Double Line Closed Blood Sampling System	. 20
TranStar® HemoDraw® Plus Triple Line Closed Blood Sampling System	. 21
LogiCal® Kids Kit™ Neonatal Closed Blood Sampling System	. 22
LogiCal® HemoDraw® Plus Single Line Closed Blood Sampling System	. 23
LogiCal® HemoDraw® Plus Double Line Closed Blood Sampling System	. 24
LogiCal® HemoDraw® Plus Triple Line Closed Blood Sampling System	. 24

Add-On Closed Blood Sampling Kits

ORDERING INFORMATION


Add-On Kids Kit™ NeoNatal Closed Blood Sampling

CodeDescriptionUnits per CaseMX413319" (48 cm) Closed Blood Sample Kit10



MX4137 1" (3 cm) Closed Blood Sample Kit

10



Add-On HemoDraw® Plus Closed Blood Sampling

CodeDescriptionUnits per CaseMX401R151HemoDraw® Plus blood sampling 152 cm/60" add on kit with one sampling site10




p. 16 Pressure Monitoring


TranStar® Kids Kit™ Neonatal Closed Blood Sampling System

ORDERING INFORMATION


TranStar® Kids Kit™ Neonatal Closed Blood Sampling System Code Description Units per Case MX9543 TranStar® Kids Kit™ 19" (48 cm) Closed Blood Sample Set (Neonatal) 10



MX9544 TranStar® Kids Kit™ 31" (79 cm) Closed Blood Sample Set (Neonatal) 10




MX9545 TranStar® Kids Kit™ 34" (86 cm) Closed Blood Sample Set (Neonatal) 10




TranStar® Kids Kit™ Neonatal Closed Blood Sampling System

ORDERING INFORMATION

TranStar®	anStar® Kids Kit™ Neonatal Closed Blood Sampling System	
Code	Description	Units per Case
MX9547	TranStar® Kids Kit™ 7" (18 cm) Closed Blood Sample Set (Neonatal)	10








p. 18 Pressure Monitoring

TranStar® HemoDraw® Plus Single Line Closed Blood Sampling System


ORDERING INFORMATION




TranStar® HemoDraw® Plus Double Line Closed Blood Sampling System

ORDERING INFORMATION


TranStar® Her	ar® HemoDraw® Plus Double Line Closed Blood Sampling System	
Code	Description	Units per Case
MX952R181	TranStar® 183 cm/72" Double Line HemoDraw® Plus blood sampling kit with one sampling site, arm mounted	10



MX952H181 TranStar® 183 cm/72" Double Line HemoDraw® Plus blood sampling kit with one sampling site 10



MX952H182 TranStar® 183 cm/72" Double Line HemoDraw® Plus blood sampling kit with two sampling sites 10




p. 20 Pressure Monitoring

TranStar® HemoDraw® Plus Triple Line Closed Blood Sampling System

ORDERING INFORMATION


TranStar® He	TranStar® HemoDraw® Plus Triple Line Closed Blood Sampling System	
Code	Description	Units per Case
MX953H182	TranStar® 183 cm/72" Triple Line HemoDraw® Plus blood sampling kit with two sampling sites	10



LogiCal® Kids Kit™ Neonatal Closed Blood Sampling System


ORDERING INFORMATION

LogiCal® Kids Kit™ Neonatal Closed Blood Sampling System Code Description Units per Case MX9644 LogiCal® Kids Kit™ 31" [76 cm] Closed Blood Sample Set (Neonatal) 10




MX9645 LogiCal® Kids Kit™ 34" (83 cm) Closed Blood Sample Set (Neonatal)

10




MX9647 LogiCal® Kids Kit™ 7" (18 cm) Closed Blood Sample Set (Neonatal)

10



MX9648 LogiCal[®] Kids Kit[™] 10" (25 cm) Closed Blood Sample Set (Neonatal)

10



p. 22 Pressure Monitoring

LogiCal® HemoDraw® Plus Single Line Closed Blood Sampling System

ORDERING INFORMATION

Code	Description	Units per Case
MX961R181	LogiCal® 183 cm/72" Single line HemoDraw® Plus blood sampling kit with one sampling site, arm mounted	10
	48" 6" 12"	
MX961H181	LogiCal® 183 cm/72" Single line HemoDraw® Plus blood sampling kit with one sampling site	10
	17" 36" 12" (C)	
ИХ961Н211	LogiCal® 213 cm/84" Single line HemoDraw® Plus blood sampling kit with one sampling site	10
	17" 48" 12" (O)	
иХ961H182	LogiCal® 183 cm/72" Single line HemoDraw® Plus blood sampling kit with two sampling sites	10
	17" 4" 32" 12"	
MX961H212	LogiCal® 213 cm/84" Single line HemoDraw® Plus blood sampling kit with two sampling sites	10
	17" 4" 43" 12"	

LogiCal® HemoDraw® Plus Double Line Closed Blood Sampling System


ORDERING INFORMATION

Code	Description	Units per Case
MX962R181	LogiCal® 183 cm/72" Double line HemoDraw® Plus blood sampling kit with one sampling site, arm mounted	10
	48" 6" 12" 48" 12"	
MX962H181	LogiCal® 183 cm/72" Double line HemoDraw® Plus blood sampling kit with one sampling site	10
	17" 36" 12" 48" 12"	
MX962H182	LogiCal® 183 cm/72" Double line HemoDraw® Plus blood sampling kit with two sampling sites	10
	17" 4" 32" 12" 48" 12"	

LogiCal® HemoDraw® Plus Triple Line Closed Blood Sampling System

ORDERING INFORMATION

LogiCal® HemoDraw® Plus Triple Line Closed Blood Sampling System			
Code	Description	Units per Case	
MX963H182	LogiCal® 183 cm/72" Triple Line HemoDraw® Plus blood sampling kit with two sampling sites	10	



p. 24 Pressure Monitoring



ACCESSORIES AND COMPONENTS





FranStar® Mounting Accessories	28
_ogiCal® Mounting Accessories	28
Blood Sampling Access Devices	29
HemoDraw® Plus Mounting Accessories	29
Flush Devices	30
Pressure Monitoring Tubing with Male/Female Luer Lock	30
Pressure Monitoring Tubing with Male/Male Luer Lock	31
Pressure Monitoring Tubing with Male/Female Luer Lock (Red)	31
Pressure Monitoring Tubing with Male/Female Luer Lock (Low Volume)	31
Non-Vented Administration Sets	32
Pressure Monitoring Components	33
Manometers	34

TranStar® Mounting Accessories

ORDERING INFORMATION

TranStar® Mounting Accessories		
Code	Description	Units per Case
MX800	TranStar® Modular Transducer Mounting Plates	20
MX240	TranStar® Pole Clamp for Mounting a Transducer Plate (Fits 3/4" to 11/2" IV Pole)	1
MX244	TranStar® Pole Clamp for Mounting a Transducer Plate (Fits 3/8" to 1" IV Pole)	1
MX800S	Armband for Patient Monitoring	10










MX244

MX800S

LogiCal® Mounting Accessories

ORDERING INFORMATION


LogiCal ® Mounting Accessories		
Code	Description	Units per Case
MX261	IV Pole Clamp with Bracket for 1 Transducer Mount Plate	1
MX262	Bracket for 2 Transducer Mounting Plates	1
MX263	Bracket for 3 Transducer Mounting Plates	1
MX264P1	Bracket for 4 Transducer Mounting Plates	1
MX260PM	Armband for Patient Mount 1 Transducer Mounting Plate	50



p. 28 Pressure Monitoring

Blood Sampling Access Devices


ORDERING INFORMATION



HemoDraw® Plus Mounting Accessories

ORDERING INFORMATION

Reservoir Bracket		
Code	Description	Units per Case
MX257C1	Reservoir Bracket for TranStar® Mounting Plate	50
MX267C1	Reservoir Bracket for LogiCal® Mounting Bracket	50



Flush Devices

ORDERING INFORMATION

Flush Devices				
Code	Description	Units per Case		
MX503A	2-4 cc/hr Flush Device with Locknut	50		
MX503GYA	2-4 cc/hr Trigger Flush	50		

Pressure Monitoring Tubing with Male/Female Luer Lock

ORDERING INFORMATION

Pressure Monitoring Tubing		
Code	Description	Units per Case
MX560	6" (15.2 cm) PVC Pressure Monitoring Line	25
MX561	12" (30.5 cm) PVC Pressure Monitoring Line	25
MX562	24" (61.0 cm) PVC Pressure Monitoring Line	25
MX563	36" (91.4 cm) PVC Pressure Monitoring Line	25
MX563L	36" (91.4 cm) PVC Pressure Monitoring Line, for use with Linden Lock Nut	25
MX564	48" (121.9 cm) PVC Pressure Monitoring Line	25
MX566	72" (182.9 cm) PVC Pressure Monitoring Line	25
MX567	84" (213.4 cm) PVC Pressure Monitoring Line	25
MX568	96" (243.8 cm) PVC Pressure Monitoring Line	25



Pressure Monitoring Tubing with Male/Male Luer Lock

ORDERING INFORMATION

Pressure Monitoring Tubing		
Code	Description	Units per Case
MX570	6" (15.2 cm) PVC Pressure Monitoring Line	25
MX571	12" (30.5 cm) PVC Pressure Monitoring Line	25
MX572	24" (61.0 cm) PVC Pressure Monitoring Line	25
MX573	36" [91.4 cm] PVC Pressure Monitoring Line	25
MX574	48" (121.9 cm) PVC Pressure Monitoring Line	25
MX576	72" (182.9 cm) PVC Pressure Monitoring Line	25
MX578	96" [243.8 cm] PVC Pressure Monitoring Line	25



p. 30 Pressure Monitoring

Pressure Monitoring Tubing with Male/Female Luer Lock (Red)

ORDERING INFORMATION

Pressure Monitoring Tubing		
Code	Description	Units per Case
MX540R	6" (15.2 cm) Guide Stripe PVC Pressure Monitoring Line	25
MX541R	12" (30.5 cm) Guide Stripe PVC Pressure Monitoring Line	25
MX542R	24" (61.0 cm) Guide Stripe PVC Pressure Monitoring Line	25
MX543R	36" (91.4 cm) Guide Stripe PVC Pressure Monitoring Line	25
MX544R	48" (121.9 cm) Guide Stripe PVC Pressure Monitoring Line	25
MX546R	72" (182.9 cm) Guide Stripe PVC Pressure Monitoring Line	25



Pressure Monitoring Tubing with Male/Female Luer Lock (Low Volume)

ORDERING INFORMATION

Pressure Monitoring Tubing				
Code	Description	Units per Case		
MX661	12" (30.5 cm) Low Volume Polyethylene Pressure Monitoring Tubing	25		
MX663	36" (91.4 cm) Low Volume Polyethylene Pressure Monitoring Tubing	25		
MX664	48" (121.9 cm) Low Volume Polyethylene Pressure Monitoring Tubing	25		
MX666	72" (182.8 cm) Low Volume Polyethylene Pressure Monitoring Tubing	25		




Non-Vented Administration Sets

ORDERING INFORMATION


Fluid Administration and Extension SetsCodeDescriptionUnits per CaseMX821SFluid Administration Set25



MX822S Bifurcated Fluid Administration Set 25



MX823S Trifurcated Fluid Administration Set 25



p. 32 Pressure Monitoring

Pressure Monitoring Components


ORDERING INFORMATION

Components		
Code	Description	Units per Case
MX43660	6" (15 cm) Low Volume PVC Pressure Monitoring Extension with Male Luer Lock and 3-Way Stopcock	25
MX43661	12" (30 cm) Low Volume PVC Pressure Monitoring Extension with Male Luer Lock and 3-Way Stopcock	25
MX43860	6" (15 cm) Low Volume PVC Pressure Monitoring Extension with Male Luer Lock and 4-Way Stopcock	25
MX801	Hydrophobic Bacterial Filter with 0.22 Micron Air Filter	50
MX805	Hydrophobic Monitoring Line Filter with 0.45 Micron Filter and Male Luer Lock	25
MX999	PA-CVP Bridge	25

Manometers

ORDERING INFORMATION

Manometers		
Code	Description	Units per Case
MX441	35 cm Manometer Set: 4-Way Stopcock, 50" Extension Set with "Y" Site 34" from Male Luer Lock, includes 3/4" and 1" Pole Clamps	10



MX447 35 cm Manometer Set with 3-Way Stopcock, Hanger Strap 20



p. 34 Pressure Monitoring

Pressure Monitoring Cables


ORDERING INFORMATION

Single, Double, Triple and Quad Channel Cables

Description

Pressure Monitoring cables to connect pressure transducers (TranStar® or LogiCal®) to patient monitors can be ordered by contacting Customer Service at 800.258.5361. Please specify the type of:

- Pressure Transducer (TranStar® or LogiCal®)
- Monitor Manufacturer Name
- Monitor Manufacturer Model Number
- Single, Double, Triple or Quad Branched Cables



For more details on our cable offerings, please refer to the cable catalog MX195050EN.

Index

CODES

MX240	28	MX666	31	MX963H182	24
MX244	28	MX800	28	MX999	33
MX257C1	29	MX800S	28	MX1821B	4
MX260PM	28	MX801	33	MX1821G	4
MX261	28	MX805	33	MX4133	16
MX262	28	MX821S	32	MX4137	16
MX263	28	MX822S	32	MX4705	4
MX264	28	MX823S	32	MX4710	4
MX264P1	28	MX920NF	29	MX4805	4
MX267C1	29	MX921NF	29	MX4810	4
MX401R151	16	MX922NF	29	MX4830	4
MX441	34	MX923NF	29	MX9501T	8
MX447	34	MX950	10	MX9502T	9
MX503A	30	MX950L	10	MX9504T	8
MX503GYA	30	MX950SCFT	10	MX9505T	8
MX540R	31	MX950SCNT	10	MX9506T	9
MX541R	31	MX950X2SC	10	MX9510T	8
MX542R	31	MX950XSC	10	MX9513T	9
MX543R	31	MX951H181	19	MX9530T	8
MX544R	31	MX951H182	19	MX9535H001C1	21
MX546R	31	MX951H211	19	MX9543	17
MX560	30	MX951H212	19	MX9544	17
MX561	30	MX951R181	19	MX9545	17
MX562	30	MX952H181	20	MX9547	18
MX563	30	MX952H182	20	MX9548	18
MX563L	30	MX952R181	20	MX9601A	11
MX564	30	MX953H182	21	MX9602A	12
MX566	30	MX960	13	MX9604A	11
MX567	30	MX960SCFA	13	MX9605A	11
MX568	30	MX960X2SC	13	MX9606A	12
MX570	30	MX960XX	13	MX9625H003C1	24
MX571	30	MX960XY	13	MX9625H004C1	24
MX572	30	MX961H181	23	MX9630A	11
MX573	30	MX961H182	23	MX9644	22
MX574	30	MX961H211	23	MX9645	22
MX576	30	MX961H212	23	MX9647	22
MX578	30	MX961R181	23	MX9648	22
MX661	31	MX962H181	24	MX43660	33
MX663	31	MX962H182	24	MX43661	33
MX664	31	MX962R181	24	MX43860	33

For more information visit our website at www.smiths-medical.com

PRODUCT(S) DESCRIBED MAY NOT BE LICENSED OR AVAILABLE FOR SALE IN CANADA AND OTHER COUNTRIES

Smiths Medical ASD, Inc. 6000 Nathan Lane North Minneapolis, MN, 55442, USA Tel: 1-614-210-7300 Toll-Free USA 1-800-258-5361 www.smiths-medical.com

MDLBCA-1072



Find your local contact information at: www.smiths-medical.com/customer-support

Smiths Medical is part of the global technology business Smiths Group plc. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. Medex, TranStar, Kids Kits, SECURE, HemoDraw, LogiCal, NovaTrans, Clear-Cuff, C-Fusor, and the Smiths Medical are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trademarks or service marks of their respective owners. Product referenced is CE-marked. ©2018 Smiths Medical. All rights reserved. SM196407EN-052018