

www.lifecorporation.com orders@lifecorporation.com

800-700-0202 | 414-272-4000

LIFE® SoftPac



LIFE-2-612

6 and 12 LPM (NORM and HIGH)

LIFE-2-025

0-25 LPM - variable flow 12 liter-flow settings

LIFE-2-6FF

6 LPM Fixed-Flow

LIFE-2-101

Replacement Cylinder M9, 250L

AED Companion O₂ Unit

Save lives by using the LIFE SoftPac Emergency Oxygen Unit to administer supplemental oxygen safely and effectively during life-threatening emergencies. This incredibly lightweight unit weighs just over 6 pounds and has a comfortable shoulder strap, making this unit great for EMTs that frequently need to carry their equipment and supplies with them. It can be used with victims of all ages, and it comes with an inhalator to use with breathing victims and includes a LIFE CPR Resuscitation mask with one-way valve for use with non-breathing victims. It features a constant reading supply gauge with simple full-to-empty symbols, always visible through clear window, knurled-knob On/Off valve prevents accidental opening in transit. In compliance with AHA regulations, this unit provides 100% inspired pure oxygen and offers a variable version of 12 different liter flow rates with click-stop settings at 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, and 25 LPMs.

- 4 1/2" x 8" x 15", 6 pounds
- 5 Year Warranty
- Inhalator for a breathing victim and resuscitator for a non-breathing victim
- Cylinder is easily refillable/replaceable
- Ready-to-use, full of oxygen
- 40 minute supply at 6 LPM

No Rx Required • Ready-to-Use • Portable • Emergency Oxygen • Refillable