

# CRUTCHES

## Aluminum, Universal and Forearm Crutches



WE CARE LIKE FAMILY™


### ALUMINUM CRUTCH FEATURES

- Lightweight anodized aluminum frame
- Push button height adjustment
- Crutch-bow ends to prevent tripping on stairs
- Comfortable, durable underarm pad and hand grip
- Rubber tip for safety

### ALUMINUM



### UNIVERSAL



### FOREARM



*Knock-down design*



### UNIVERSAL CRUTCH FEATURES

- 3 in 1—adjusts from youth to tall sizes
- Handgrip adjusts easily with euro style clip
- Comfortable underarm pad and hand grip
- Retail friendly packaging and knock-down design

### FOREARM CRUTCH FEATURES

- Ergonomically contoured arm cuffs
- Adjustable leg and forearm sections
- Vinyl hand grips

Reorder No.	Product Name	Size	Height Range	Quantity	HCPCS
10100-8	Aluminum Crutches	Child	4'0" - 4'6"	8/Cs or 1/Cs	E0114
10101-8	Aluminum Crutches	Youth	4'6"-5'2"	8/Cs or 1/Cs	E0114
10102-8	Aluminum Crutches	Medium/Adult	5'2"-5'10"	8/Cs or 1/Cs	E0114
10103-8	Aluminum Crutches	Tall/Adult	5'10"-6'6"	8/Cs or 1/Cs	E0114
10104	Universal Crutches	Youth-Tall	4'7"-6'7"	1/Cs	E0114 or E0116
10110	Forearm Crutches	Youth	4'2" - 5'2"	6/Cs or 1/Cs	E0114 or E0116
10111	Forearm Crutches	Adult	5'0" - 6'2"	6/Cs or 1/Cs	E0114 or E0116
10112	Forearm Crutches	Tall	5'10" - 6'6"	6/Cs or 1/Cs	E0114 or E0116



Find Us On:



**dynarex®** is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR  
Call us at: (888) 396-2739 • [www.dynarex.com](http://www.dynarex.com)