



BARIATRIC SHOWER CHAIR WITH BACK

REORDER NO. 10321 – 1 PER BOX

REORDER NO. 10321-4 – 4 BOXES, 1 PER BOX

BARIATRIC SHOWER CHAIR WITHOUT BACK

REORDER NO. 10325 – 1 PER BOX

REORDER NO. 10325-4 – 4 BOXES, 1 PER BOX

BARIATRIC SHOWER CHAIRS

500 LB. WEIGHT CAPACITY

- Constructed with a corrosion-proof aluminum frame and blow-molded plastic seat
- Cross-brace design provides additional support
- Seat features drainage holes and integrated hand grips for added security
- Angled legs with rubber tips for added security
- Lightweight for easy movement
- Tool-free assembly and height adjustment in 1-inch increments
- 10321: Blow-molded plastic back is removable for added convenience
- Comes in retail-friendly packaging
- Limited Lifetime Warranty
- HCPCS: E0240 (10321), E0245 (10325)



10321



10325