

# DYNAREST™ FOAM MATTRESSES

## Premium Visco-Gel Infused Foam Mattress





WE CARE LIKE FAMILY™

### FEATURES AND BENEFITS

- Constructed from three layers of variable density foam
- Top layer features visco-gel infused memory foam for improved pressure redistribution
- Middle layer features perpendicular air channels for improved air flow
- Firm bottom layer provides maximum support
- Removable, fluid resistant cover is easy to maintain, reduces perspiration, and provides low shear friction
- Meets Flame Retardant CFR 1632 and CFR 1633 standards
- Weight capacity 450 lb.
- **HCPCS: E0184**

### PERIMETER OPTIONS AVAILABLE



Reorder No.	Product Name	Quantity
10425	DynaRest Premium Visco-Gel Infused Foam Mattress, 36" x 76"	1/Bx
10426	DynaRest Premium Visco-Gel Infused Foam Mattress, 36" x 80"	1/Bx
10427	DynaRest Premium Visco-Gel Infused Foam Mattress, 36" x 84"	1/Bx
10428	DynaRest Premium Visco-Gel Infused Foam Mattress w/Perimeter, 36" x 80"	1/Bx



Find Us On:



**dynarex®** is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR  
Call us at: (888) 396-2739 • [www.dynarex.com](http://www.dynarex.com)