



## REORDER NUMBERS:

10445 – 42"×80"

10446 – 48"×80"

10439 – 54"×80"

## BARIATRIC HD AIRFLOAT AIR MATTRESS WITH PUMP

- Designed to help prevent and relieve pressure wounds
- Equipped with sixteen 10" pressure cycling air cells, and four static cells at the head section
- Features four inflation modes: alternating cycle, static mode, seat function, and max firm
- Alternating pressure helps to redistribute weight for improved circulation
- Static mode provides steady, consistent pressure throughout all air cells
- Seat function adjusts the alternating air cells for ideal support in a seated position
- Max firm feature fills the air cells to the maximum firmness for safe transfers
- Digital button pump features 10, 15, 20 and 25-minute cycle times
- Audible alarm sounds to alert in case of power failure or malfunction
- Limited Warranty: 24 Months
- HCPCS: E0277



CPR VALVE ALLOWS FOR RAPID DEFLATION



## WEIGHT CAPACITY

- **850 lb.** (10445-10446)
- **1000 lb.** (10439)