

# DYNAREST™ AIR MATTRESSES

## Airfloat 100 Air Mattress with Pump



WE CARE LIKE FAMILY™

### FEATURES & BENEFITS

- Alternating pressure mattress replacement system
- Twenty 8-inch deep air cells with eighteen laser holes provide low air loss
- Pressure control dial with weight indicators adjust pressure levels
- CPR valve allows for rapid deflation
- 2 modes: Alternating cycle and Static
- Equipped with low pressure indicator light
- **HCPCS: E0277**



Reorder No.	Product Name	Dimensions	Cell Height	Cell on Cell	Cover Type	Max. Capacity
10441	DynaRest Airfloat 100 36" Air Mattress with Pump	36"x80"	8"	–	Nylon	350 lb.
10448	DynaRest Airfloat 100 36" Long Air Mattress with Pump	36"x84"	8"	–	Nylon	400 lb.
10449	DynaRest Airfloat 100F 30" Air Mattress with Foam Base and Pump	30"x80"	5" air, 3" foam	foam base	Nylon	400 lb.



Find Us On:



**dynarex®** is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • [www.dynarex.com](http://www.dynarex.com)