

HYDRAULIC PATIENT LIFT

Transfer Patients Who Have Limited Mobility



WE CARE LIKE FAMILY™

FEATURES AND BENEFITS

- High performance hydraulics gradually raise and lower patients
- Maximum weight capacity 450 lb.
- Durable swivel bar accommodates 2–6 sling straps
- Sling hook-up height adjusts between 16"–60"
- Compatible with standard slings—with and without commode opening, chains, and straps
- Packaged with Chains (10753)
- 5-year structural steel warranty, 1-year mechanical parts/components warranty
- Tool-free assembly



Reorder No.	Product Name	Quantity
10704	Hydraulic Patient Lift with Chains	1/Cs

COMPARE TO:

Painted Hydraulic Lift™ by Invacare

Hydraulic Deluxe Silver Vein Patient Lift™ by Drive DeVilBiss Healthcare

ALSO AVAILABLE:

Lift Scale (10705) Standard Canvas (10750CN-10751CN) and Mesh Slings (10750MS-10751MS)
(with and without commode opening)

Reorder No.	Product Name	Quantity
10705	Lift Scale	1/Bx
10750MS	Standard Mesh Sling	1/Bx
10750CN	Standard Canvas Sling	1/Bx
10751MS	Standard Mesh Sling with Commode Opening	1/Bx
10751CN	Standard Canvas Sling with Commode Opening	1/Bx
10753	Replacement Chains for Standard Slings	1/Bx
10754	Replacement Straps for Standard Slings	1/Bx



Find Us On:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • www.dynarex.com