



10740XX



10741XX



10742XX



10743XX

LIFT SLINGS

FOR USE WITH MOST ELECTRIC PATIENT LIFTS

- Designed to assist in safe patient transfers from bed, wheelchair, or floor
- 10742XX: Transports patients for toileting purposes
- Manufactured from strong, durable fabric for reliable support
- 10740XX & 10742XX: Features breathable mesh material for optimal comfort
- Designed to fit most electric patient lifts for optimal versatility
- Features convenient back strap for easier positioning
- Integrated color-coded hooks help assure correct patient positioning
- Color-coded perimeter bindings allow for easy size-identification
- Ideal for use in a variety of long-term care settings
- Six Month Limited Warranty
- HCPCS: E0621



10740XXXL

WEIGHT CAPACITY

- **600 lb.** (10740XX, 10741XX, 10742XX, 10743XX)
- **1000 lb.** (10740XXXL)

BARIATRIC FULL BODY SLING

10740XX	MESH	XX-LARGE
10740XXXL	MESH	XXX-LARGE
10741XX	SOLID FABRIC	XX-LARGE

BARIATRIC FULL BODY SLING WITH COMMODE OPENING

10742XX	MESH	XX-LARGE
---------	------	----------

BARIATRIC DIVIDED LEG SLING

10743XX	SOLID FABRIC	XX-LARGE
---------	--------------	----------

