

LONGTERM CARE BEDS

D500 Series Ultra Low Beds



WE CARE LIKE FAMILY™

FEATURES AND BENEFITS

- Constructed from welded steel
- Features true verticle lift
- Deck to floor range 2.5"—26"
- Split frame design along with a transport bar transforms bed into a convenient to move unit
- Integrated wooden head and foot boards
- Ultra-quiet DC motors with 5 functions including Trendelenburg, and Reverse Trendelenburg
- Equipped with four 3" casters



Warranty: 15 years on structural steel, 3 years on electrical & mechanical, 1 year on all other parts and components from original date of purchase.

Transport Mode



*Low Position,
2" off Floor*



Series	Height		Deck	Extended Heel Lift	Weight Capacity
	Min	Max			
D500	2.5"	26"	Grid	✓	400 lb.



Mahogany Amber Cherry Light Oak



Find Us On:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR
Call us at: (888) 396-2739 • www.dynarex.com