



CLOSED SUCTION CATHETERS TRACHEAL

DOUBLE SWIVEL ELBOW • LENGTH: 30.5 CM/12 IN



Dynarex Resp-O₂™ Closed Suction Endotracheal Catheters are disposable respiratory devices designed to safely remove secretions from the respiratory tract of ventilator-dependent patients. The device allows ventilation to continue during suctioning, helps prevent cross-contamination, and protects caregivers from airborne pathogens.

- Protects patient from contact with outside pathogens and prevents contamination of the surrounding environment
- Smooth, flexible tubing minimizes patient discomfort
- Pressure wiper seal helps maintain PEEP
- Flex Tube provides adjustment and flexibility to the ventilator circuit
- Designed for 24-hr, single-patient use
- Does not contain DEHP

10 Fr

Double Swivel Elbow

REORDER NO. 31010

ADULT/PEDIATRIC

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

12 Fr

Double Swivel Elbow

REORDER NO. 31011

ADULT/PEDIATRIC

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

14 Fr

Double Swivel Elbow

REORDER NO. 31012

ADULT

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

16 Fr

Double Swivel Elbow

REORDER NO. 31013

ADULT

ACCESSORIES

15 mm x 22 mm Flex Tube
Change of Day

PACKAGE QUANTITY
20/Case

