


CO₂ MASKS & ACCESSORIES

Dynarex Resp-O₂™ Oxygen CO₂ Masks & Accessories are designed to help monitor patients' exhaled carbon dioxide levels while maintaining the continuous administration of oxygen. The masks feature an elongated shape for added patient comfort and swivel design for optimal flexibility and range of motion.

- Designed to monitor CO₂ levels from the nose and mouth
- Soft, elongated mask provides a secure fit
- Transparent materials allow for clear visualization
- Includes 7' of crush-resistant, star lumen oxygen supply tubing
- 33080 & 33081 feature CO₂ sampling line and O₂ tubing of equal length
- Not made with natural rubber latex or DEHP



**Equal Length
CO₂ Sampling
Line & O₂
Tubing**

33080 & 33081



Oxygen CO₂
Elongated Mask

**7' (2.1 M) CO₂ SAMPLING LINE
WITH FILTER
7' (2.1 M) OXYGEN TUBE
MALE/MALE CONNECTOR**

REORDER NO. 33080

ADULT

PACKAGE QUANTITY
25/Case

Oxygen CO₂
Elongated Mask

**7' (2.1 M) CO₂ SAMPLING LINE
WITH FILTER
7' (2.1 M) OXYGEN TUBE
MALE/MALE CONNECTOR**

REORDER NO. 33081

PEDIATRIC

PACKAGE QUANTITY
25/Case

Oxygen CO₂ 7' (2.1m)
Sampling Line

FEMALE/FEMALE CONNECTOR

REORDER NO. 33082

ALL AGES

PACKAGE QUANTITY
25/Case

Oxygen CO₂ 7' (2.1m)
Sampling Line

MALE/MALE CONNECTOR

REORDER NO. 33083

ALL AGES

PACKAGE QUANTITY
25/Case

Oxygen CO₂ 7' (2.1m)
Sampling Line

MALE/FEMALE CONNECTOR

REORDER NO. 33084

ALL AGES

PACKAGE QUANTITY
25/Case

Oxygen CO₂ 7' (2.1m)
Sampling Line

**MALE/MALE CONNECTOR
WITH FILTER**

REORDER NO. 33085

ALL AGES

PACKAGE QUANTITY
25/Case

Oxygen CO₂ 7' (2.1m)
Sampling Line

**MALE/FEMALE CONNECTOR
WITH FILTER**

REORDER NO. 33086

ALL AGES

PACKAGE QUANTITY
25/Case

Oxygen CO₂
Pressure Line Filter

REORDER NO. 33087

ALL AGES

PACKAGE QUANTITY
50/Case

