



NASAL OXYGEN CANNULAS Flared Tip

Dynarex Resp-O₂™ Nasal Oxygen Cannulas with Flared Tips are designed to efficiently deliver supplemental oxygen to respiratory patients. The cannulas are available with either Standard or Universal connectors and feature light green oxygen supply tubing, designed to help improve tube visibility.

- Flared-tip nasal prongs help ensure a secure fit
- Crush-resistant, star lumen oxygen supply tubing
- Individually wrapped
- Not made with natural rubber latex or DEHP



Star Lumen



Universal Connector



Standard Connector



7' (2.1 m) Crush Resistant Tubing

STANDARD CONNECTOR

REORDER NO. 33110

ADULT

PACKAGE QUANTITY
50/Case

7' (2.1 m) Crush Resistant Tubing

UNIVERSAL CONNECTOR

REORDER NO. 33111

ADULT

PACKAGE QUANTITY
50/Case

14' (4.3 m) Crush Resistant Tubing

STANDARD CONNECTOR

REORDER NO. 33112

ADULT

PACKAGE QUANTITY
50/Case

21' (6.4 m) Crush Resistant Tubing

STANDARD CONNECTOR

REORDER NO. 33113

ADULT

PACKAGE QUANTITY
25/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210426