



AEROSOL MASKS

Dynarex Resp-O₂™ Aerosol Masks are designed for use in aerosol therapy to help administer inhalable treatments and aerosol medications to respiratory patients.

- Transparent vinyl mask features smooth edges
- Mask features an adjustable nose clip and elastic strap
- Can be used with most handheld nebulizers
- Connector accepts 22 mm aerosol corrugated tubing and nebulizer cup
- Not made with natural rubber latex or DEHP



ELONGATED



FACE TENT MASK



STANDARD



Aerosol Standard Mask

REORDER NO. 34000

ADULT

PACKAGE QUANTITY
50/Case

Aerosol Standard Mask

REORDER NO. 34001

PEDIATRIC

PACKAGE QUANTITY
50/Case

Aerosol Standard Mask

REORDER NO. 34002

INFANT

PACKAGE QUANTITY
50/Case

Aerosol Elongated Mask

REORDER NO. 34003

ADULT

PACKAGE QUANTITY
50/Case

Aerosol Elongated Mask

REORDER NO. 34004

PEDIATRIC

PACKAGE QUANTITY
50/Case

Aerosol Face Tent Mask

REORDER NO. 34005

ADULT

PACKAGE QUANTITY
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

©2021 Dynarex Corporation

Contact your local sales representative or call to find your nearest distributor.



R210427