



OXYGEN TANK CARTS AND HOLDERS

Dynarex Resp-O₂™ Oxygen Tank Carts are designed to safely and reliably carry and support D- and E-sized oxygen cylinders. The sturdy frame is constructed from durable carbon steel with an 80-lb capacity. The cylinder cart is available in single-tank and dual-tank varieties and features a height-adjustable handle for easy customization.

- Durable yet lightweight cart provides optimal portability
- 36620 & 36621: Silver-vein finish is easy to maintain
- 36622: Chrome plating helps to resist corrosion
- Large, composite wheels allow for smoother maneuverability
- 80-lb. weight capacity
- Not made with natural rubber latex or DEHP



Oxygen Tank Cart

REORDER NO. 36620

PACKAGE QUANTITY
1/Case



Dual Oxygen Tank Cart

REORDER NO. 36621

PACKAGE QUANTITY
1/Case



Oxygen Tank Holder Wheelchair Mount

REORDER NO. 36622

PACKAGE QUANTITY
1/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210823