



WE CARE LIKE FAMILY™

DENTAL PRODUCTS



5 OZ. DRINKING CUPS

Reorder No. 4236 – White

Reorder No. 4237 – Blue

Reorder No. 4238 – Mint Green

Reorder No. 4239 – Mauve

Reorder No. 4240 – Lavender

- Constructed from crush-resistant polypropylene for optimal sturdiness
- Rolled lip allows for easy separation
- Ribbed center provides added grip and enhanced durability
- Convenient 5 oz. size for easy dispensing
- 50 per box, 20 boxes per case



DENTAL BIBS

Reorder No. 4491 – White

Reorder No. 4492 – Blue


Reorder No. 4493 – Mint Green

Reorder No. 4494 – Mauve

Reorder No. 4495 – Lavender

Reorder No. 4496 – Black

- Standard size approximately 13" × 18" to provide ample coverage
- 2 ply-paper provides optimal absorbency
- Water-resistant polyethylene backing helps to prevent soak-through
- 500 per case



www.dynarex.com • 888-396-2739

Dynarex Corporation • 10 Glenshaw Street • Orangeburg, NY 10962

R211104



© 2022 Dynarex Corporation