



WE CARE LIKE FAMILY™

PLASTIC THUMB FORCEPS/ TWEEZERS 4.25"



**Ideal for grasping items from surfaces.
Used for medical and first aid purposes.**

- Tapered ends for precise pick-ups
- Sturdy and disposable
- Unique construction maintains tip alignment
- Straight serrated tip
- Non-locking
- Anti-slip thumb and finger ridges for extra control
- Not made with natural rubber latex



Reorder No. 4917: 3,000 per case (bulk packed)

Reorder No. 4918: 1,000 per case (Individually wrapped)

Also available:

Reorder Nos. 3102-3134: Stretch Gauze Bandages

Reorder Nos. 3651-3657: Waterproof Adhesive Tape in a Plastic Spool

www.dynarex.com • 888-396-2739

Dynarex Corporation • 10 Glenshaw Street • Orangeburg, NY 10962

R210820



© 2021 Dynarex Corporation