

The Brand to Know!™



Quality Medical Products Since 1967



# Urethral Catheterization Tray

For convenient and economical catheterization procedures



4925



4928

Reorder No.	Description	Quantity
4925	Urethral Catheterization Tray	20/Cs

Reorder No.	Description	Quantity
4928	Urethral Catheterization Tray	20/Cs

### Includes:

- |   |   |
|---|---|
| 1 Graduated Collection Basin            | 3 Povidone-Iodine Swabsticks                |
| 1 Specimen Container With Lid And Label | 1 Pair Cuffed Powder-Free Vinyl Exam Gloves |
| 1 14Fr. Red Rubber Catheter             | 3 Underpads                                 |

- Components packed within tray to assure components are available in proper sequence during procedure
- Sterile



Find us on:



### Includes:

- |                                       |   |
|---------------------------------------|---|
| 1 14 Fr. (Male) Intermittent Catheter | 1 Waterproof Drape (Blue, 45x42 cm)         |
| 1 1,000cc Outer Basin Tray            | 1 Lubricating Jelly Packet (3g)             |
| 1 4oz. Specimen Cup                   | 3 Povidone-Iodine Swabsticks                |
| 1 Fenestrated Drape (White, 45x42 cm) | 1 Pair Cuffed Powder-Free Vinyl Exam Gloves |
| 1 Patient Label                       |   |

- Straight tip
- Medical grade PVC; not made with natural rubber latex
- Offset polished eyelets to maximize urine flow
- Sterile

See our website for Urinary Drainage Bags (#4270-4271)



For more information contact your local Dynarex® distributor sales representative.

Or call 888-396-2739 to find your nearest distributor.