





## Sterile Intermittent Catheters



### Trust in Dynarex for Patient Care

Dynarex Intermittent Straight Catheters effectively manage bladder drainage in cases of urinary retention, bladder obstruction, paralysis or loss of sensation caused by spinal cord injury. They are available in male and female styles, and are for single use only. Regular draining of the bladder will protect the patient's kidneys from infection and will lower the risk of distending the bladder.



- **Straight tip**
- **Latex-free**
- **Color-coded for size visualization**
- **Medical grade PVC**
- **Offset polished eyelets to maximize urine flow**
- **Individually sealed in easy-open, peel-strip sterile packages**



Reorder No.	Intermittent Straight Catheters	Quantity
4981	Male – 8 French, Blue	50/cs
4982	Male – 10 French, Black	50/cs
4983	Male – 12 French, White	50/cs
4984	Male – 14 French, Green	50/cs
4985	Male – 16 French, Orange	50/cs
4986	Male – 18 French, Red	50/cs
4991	Female – 8 French, Blue	50/cs
4992	Female – 10 French, Black	50/cs
4993	Female – 12 French, White	50/cs
4994	Female – 14 French, Green	50/cs
4995	Female – 16 French, Orange	50/cs

Also see our catalog for Urinary Drainage Bags (#4270-4271) and Urinary Leg Bags (#4281-4282)



Depending on market demand and Dynarex stock, we will be happy to address any orders for items/sizes currently not shown.

© 2012 Dynarex Corporation