



WE CARE LIKE FAMILY™

STERILE NITRILE EXAM GLOVES



SINGLE

Reorder No. 6522 – Small

Reorder No. 6523 – Medium

Reorder No. 6524 – Large

PAIRS

Reorder No. 6528 – Small

Reorder No. 6529 – Medium

Reorder No. 6530 – Large

- Sterile gloves packaged in individual peel open packs
- Textured palms and beaded cuffs
- 2 mil minimum thickness
- Ambidextrous and powder-free
- Not made with natural rubber latex
- Single: 100 gloves per box, 8 boxes per case
Pairs: 50 pairs of gloves per box, 8 boxes per case

