



WE CARE LIKE FAMILY



## Insulin Syringes

**Clear, bold markings assure dosage accuracy**

- The safety syringes are designed with a safety feature which retracts the needle back into the barrel and breaks the plunger to prevent needle stick injury
- Hands and fingers remain away from the needle during injection and activation of safety mechanism
- Needle is secured into barrel of syringe, cannot be reused
- Available in .5 cc and 1 cc
- Permanent attached needle
- Not made with natural rubber latex
- Non-Toxic
- Sterile
- Non-Pyrogenic

6901-6905



6940-6959



## Pen Needles

- Used in conjunction with injection pens to deliver injectable medications into the body
- Come in a variety of needle lengths and diameters used by both health professionals and patients
- Non-Toxic, Non-Pyrogenic
- Not made with natural rubber latex
- Sterile
- Available in 30 G-5 mm, 31 G-8 mm, 30 G-8 mm, 31 G-5 mm, and 32 G-4 mm

7024-7028



HCPCS Codes: 6901-6905 A4206  
HCPCS Codes: 6940-6959 A4206  
HCPCS Codes: 7026-7028 A4215



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)

© 2016 Dynarex Corporation

# dynarex INSULIN SYRINGES AND PEN NEEDLES



## Safety Insulin Syringes With Needle

Reorder No.	Description	Quantity
6901	Safety Insulin Syringe With Needle, 0.5 cc, 29 G x 1/2"	5/100/Cs
6902	Safety Insulin Syringe With Needle, 0.5 cc, 30 G x 5/16"	5/100/Cs
6903	Safety Insulin Syringe With Needle, 1 cc, 27 G x 1/2"	5/100/Cs
6904	Safety Insulin Syringe With Needle, 1 cc, 29 G x 1/2"	5/100/Cs
6905	Safety Insulin Syringe With Needle, 1 cc, 30 G x 1/2"	5/100/Cs

## Pen Needles

Reorder No.	Description	Quantity
7024	<b>NEW!</b> Pen Needles, 30 G x 5 mm	12/100/Cs
7025	<b>NEW!</b> Pen Needles, 31 G x 8 mm	12/100/Cs
7026	Pen Needles, 30 G x 8 mm	12/100/Cs
7027	Pen Needles, 31 G x 5 mm	12/100/Cs
7028	Pen Needles, 32 G x 4 mm	12/100/Cs

Compatible with these insulin pens:

- Autopen®
- BD pen
- FlexPen®
- NovoPen® / NovoLet® / InnoLet®
- HumaPen Ergo®
- HumaPen® Luxura™
- HumaPen® Memoir™
- KwikPen™
- Solostar®
- SymlinPen®
- BerliPen®
- Omnican Pen®
- OptiClik™
- OptiPen® / OptiSet®



## Insulin Syringes With Needle (Non-Safety)

Reorder No.	Description	Quantity
6940	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 27 G x 1/2", Individually Packaged	5/100/Cs
6941	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 28 G x 1/2", Individually Packaged	5/100/Cs
6942	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 29 G x 1/2", Individually Packaged	5/100/Cs
6943	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 30 G x 5/16", Individually Packaged	5/100/Cs
6944	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 31 G x 5/16", Individually Packaged	5/100/Cs
6945	Insulin Syringe With Needle, Non-Safety, 1 cc, 27 G x 1/2", Individually Packaged	5/100/Cs
6946	Insulin Syringe With Needle, Non-Safety, 1 cc, 28 G x 1/2", Individually Packaged	5/100/Cs
6947	Insulin Syringe With Needle, Non-Safety, 1 cc, 29 G x 1/2", Individually Packaged	5/100/Cs
6948	Insulin Syringe With Needle, Non-Safety, 1 cc, 30 G x 1/2", Individually Packaged	5/100/Cs
6949	Insulin Syringe With Needle, Non-Safety, 1 cc, 31 G x 5/16", Individually Packaged	5/100/Cs

Reorder No.	Description	Quantity
6950	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 27 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6951	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 28 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6952	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 29 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6953	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 30 G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs
6954	Insulin Syringe With Needle, Non-Safety, 0.5 cc, 31 G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs
6955	Insulin Syringe With Needle, Non-Safety, 1 cc, 27 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6956	Insulin Syringe With Needle, Non-Safety, 1 cc, 28 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6957	Insulin Syringe With Needle, Non-Safety, 1 cc, 29 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6958	Insulin Syringe With Needle, Non-Safety, 1 cc, 30 G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6959	Insulin Syringe With Needle, Non-Safety, 1 cc, 31 G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs