



WE CARE LIKE FAMILY™



NEW

SecureSafe™ Safety Syringes

Secure • Safe • Easy To Use

- Insulin syringes
- Tuberculin Safety Syringes
- Allergy Safety Syringe Tray



New SecureSafe™ Safety Mechanism Protects the end-user from accidental needle stick injuries

It's Easy, It's Safe, and Simple to Use!

SecureSafe™ Insulin Syringes, Tuberculin Safety Syringes and Allergy Safety Syringe Tray

The SecureSafe™ engineered safety mechanism features a cylinder safety cover to conceal the needle with a quick pull and twist movement to lock.

1 Open

2 Inject

3 Pull Up

4 Twist

See Our Website For Video Demos



See Reverse Side



Find us on:



dynarex® is your one-stop medical products brand since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

NEW COMPLETE LINE OF SAFETY SYRINGES

SecureSafe™ Safety Syringes by Dynarex

Insulin Syringes • Tuberculin Syringes • Allergy Safety Syringe Tray

**See Our Website
For Video Demos**



Reorder No.	Description	Quantity
8901	SecureSafe Safety Insulin Syringe - 0.5cc, 29G, 1/2" needle	5/100/Cs
8904	SecureSafe Safety Insulin Syringe - 1cc, 29G, 1/2" needle	5/100/Cs

Reorder No.	Description	Quantity
8937	SecureSafe Tuberculin Safety Syringe - 1cc, 25G, 5/8" needle	5/100/Cs
8938	SecureSafe Tuberculin Safety Syringe - 1cc, 26G, 3/8" needle	5/100/Cs
8939	SecureSafe Tuberculin Safety Syringe - 1cc, 27G, 1/2" needle	5/100/Cs

Reorder No.	Description	Quantity
8936	SecureSafe Allergy Safety Syringe Tray - 27G x 1/2" 1cc - 25/Tray	40/Cs

* SecureSafe Insulin, Tuberculin and Allergy syringes are comparable to BD Safety-Lok™ Syringe and Covidien Monoject™ Insulin Safety Syringe. BD Safety-Lok™ Syringe is a registered trademark of Becton, Dickinson and Company. Covidien Monoject™ Insulin Safety Syringe is a registered trademark of Medline Industries, Inc.

Now With Our New SecureSafe™ Safety Mechanics

See Front Side For Details

See our website for all our medical products! www.dynarex.com