

Surgical Sutures



WE CARE LIKE FAMILY™

Used to hold body tissues together after an injury or surgery. Application generally involves using a needle with an attached length of thread.

- Available in various thread thickness and needle dimensions to accommodate a variety of procedures
- Not made with natural rubber latex • Sterile
- All boxes individually shrink wrapped

Reorder No.	Synthetic Nylon Sutures Description	Needle Parameters and Type	Dynarex Needle Symbol	Dynarex Needle Size	Quantity
9100	Black Monofilament, Non-Absorbable, 4-0, Length 18"	PC31 - Precision REVERSE CUT 3/8 circle	▽	19 mm	12/Box
9101	Black Monofilament, Non-Absorbable, 6-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
9102	Black Monofilament, Non-Absorbable, 5-0, Length 18"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box
9103	Black Monofilament, Non-Absorbable, 4-0, Length 18"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box
9104	Black Monofilament, Non-Absorbable, 3-0, Length 18"	C7 - REVERSE CUT 3/8 circle	▼	24 mm	12/Box
9105	Black Monofilament, Non-Absorbable, 6-0, Length 18"	C1 - Precision REVERSE CUT 3/8 circle	▽	11 mm	12/Box
9106	Black Monofilament, Non-Absorbable, 5-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
9107	Black Monofilament, Non-Absorbable, 4-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
Reorder No.	Synthetic Polypropylene Sutures Description	Needle Parameters and Type	Dynarex Needle Symbol	Dynarex Needle Size	Quantity
9117	Blue Monofilament, Non-Absorbable, 4-0, Length 18"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box
9118	Blue Monofilament, Non-Absorbable, 6-0, Length 18"	C1 - Precision REVERSE CUT 3/8 circle	▽	11 mm	12/Box
9119	Blue Monofilament, Non-Absorbable, 5-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
9120	Blue Monofilament, Non-Absorbable, 4-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
Reorder No.	Synthetic Braided Polyglycolic Acid (PGA) Sutures Description	Needle Parameters and Type	Dynarex Needle Symbol	Dynarex Needle Size	Quantity
9130	Violet, Braided, Coated Synthetic Absorbable, 4-0, Length 30"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box
9131	Violet, Braided, Coated Synthetic Absorbable, 5-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
9132	Violet, Braided, Coated Synthetic Absorbable, 4-0, Length 18"	C3 - Precision REVERSE CUT 3/8 circle	▽	13 mm	12/Box
Reorder No.	Natural Material Braided Silk Sutures Description	Needle Parameters and Type	Dynarex Needle Symbol	Dynarex Needle Size	Quantity
9142	Black, Non-Absorbable, 4-0, Length 18"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box
9143	Black, Non-Absorbable, 3-0, Length 18"	C6 - REVERSE CUT 3/8 circle	▼	19 mm	12/Box



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com



WE CARE LIKE FAMILY™

Surgical Sutures



Synthetic Nylon Sutures

- A non-absorbable suture offering high tensile strength
- Monofilament for a smooth passage through tissue




Synthetic Polypropylene Sutures

- A non-absorbable suture offering high tensile strength
- Monofilament for a smooth passage through tissue



Synthetic Braided Polyglycolic Acid (PGA) Sutures

- An absorbable, violet suture made from polyglycolic acid
- Braided for strength and easy handling and coated to allow for easy passage through tissue



Natural Material Braided Black Silk Sutures

- A non-absorbable suture made from natural silk that is degummed to remove extraneous materials
- Braided for strength and easy handling



Needle Compatibility Chart

Needle Parameters	Dynarex	CP Medical	Ethicon	Covidien	LOOK
Precision REVERSE CUT, 3/8 circle, 19 mm	PC31	NPS-2	PS-2	P-12	PC31
Precision REVERSE CUT, 3/8 circle, 13 mm	C3	NP-3	P-3	P-13	C3
REVERSE CUT, 3/8 circle, 19 mm	C6	NFS-2	FS-2	C-13	C6
REVERSE CUT, 3/8 circle, 24 mm	C7	NFS-1	FS-1	C-14	C7
Precision REVERSE CUT, 3/8 circle, 11 mm	C1	NP-1	P-1	P-10	C1



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com