


Advanced Wound Care Dressings



WE CARE LIKE FAMILY™

CuraFoam™ Foam Dressing




- Maintains a moist environment, while preventing bacterial contamination
- Used to manage chronic and acute wounds—from ulcers to donor sites

COMPARE TO: HydraFoam, Medline OptiFoam & Kendall Hydrasorb

Reorder No.	Description	Quantity	HCPCS
3011	2"×2"	12/10/Cs	A6209
3012	4"×4¼"	12/10/Cs	A6210
3013	6"×6"	12/10/Cs	A6210

DynaDerm™ Hydrocolloid Dressing



- Maintains a moist wound environment and can be worn for several days
- Used on non-infected wounds of partial and full thickness

COMPARE TO: DermaRite, DermaFilm, Kendall Ultec, ConvaTec Duoderm, Medline ExuDerm

Reorder No.	Description	Quantity	HCPCS
3014	2"×2" – Thin	12/20/Cs	A6234
3015	4"×4" – Thin	12/10/Cs	A6234
3016	6"×6" – Thin	12/5/Cs	A6238
3017	6"×7" – Sacral	12/5/Cs	A6238
3018	4"×4" – Extra Thin	12/10/Cs	A6234
3019	6"×6" – Extra Thin	12/5/Cs	A6235
3020	4"×4" – Standard	12/5/Cs	A6209

DynaLevin™ Waterproof Adhesive Bordered Foam



- Maintains a moist wound environment with medium absorbency.
- Used to manage chronic and acute ulcers

COMPARE TO: DermaRite, DermaLevin, Smith & Nephew Allevyn Adhesive

Reorder No.	Description	Quantity	HCPCS
3021	4"×4"	12/10/Cs	A6212
3022	6"×6"	12/10/Cs	A6212

FoamFlex™ Non-Adhesive Waterproof Foam




- A semi-permeable film prevents bacterial contamination
- Used for ulcers and low to moderate exuding wounds

COMPARE TO: DermaRite, DermaFoam, Smith & Nephew Allevyn Foam

Reorder No.	Description	Quantity	HCPCS
3023	2"×2"	12/10/Cs	A6209
3024	4"×4"	12/10/Cs	A6209
3025	6"×6"	12/10/Cs	A6210

DynaGinate™ Calcium Alginate Dressing




- Used as a topical wound dressing for moderate to high exuding wounds
- Absorbs up to 17 times its own weight

COMPARE TO: DermaRite, DermaGinate

Reorder No.	Description	Quantity	HCPCS
3026	2"×2"	12/10/Cs	A6196
3027	4¼"×4¼"	12/10/Cs	A6197
3028	4"×8"	12/5/Cs	A6197
3029	2g × 30 cm	12/5/Cs	A6199

DynaGinate™ AG Silver Calcium Alginate Dressing



- Provides 3 days of continuous antibacterial protection with silver
- Used on high risk and infected wounds

COMPARE TO: DermaRite, DermaGinate AG, ConvaTec Aquacel AG, J&J Silvercel

Reorder No.	DynaGinate™ AG	Quantity	HCPCS
3030	2"×2"	10/10/Cs	A6196
3031	4"×5"	10/10/Cs	A6197
3032	4"×8"	10/5/Cs	A6197
3033	¾"×12" – Rope Dressing	10/5/Cs	A6199

DynaGuard™ Waterproof Composite Dressing



- Multi-layer dressing with a low adherent layer for protection
- Used on ulcers, chronic and acute wounds, and donor sites

COMPARE TO: DermaRite, DermaDress, Smith & Nephew CovSite Plus

Reorder No.	DynaGuard™	Quantity	HCPCS
3035	4"×4"	12/10/Cs	A6203
3036	6"×6"	12/10/Cs	A6203

DynaFoam™ Waterproof Bordered Foam




- Multi-layer dressing with a low adherent layer for protection
- Used on ulcers, chronic and acute wounds, and donor sites

COMPARE TO: DermaRite Bordered Foam, Medline Optifoam

Reorder No.	DynaFoam™	Quantity	HCPCS
3037	4"×4"	12/10/Cs	A6212
3038	6"×6"	12/10/Cs	A6213

Petroleum Gauze Dressing



- Non-adhering and maintains a moist wound environment
- Cotton gauze impregnated with petroleum conforms to wound
- Used for uninfected, non-draining wounds including tube sites, skin grafts, burns, and ulcers

COMPARE TO: Curad Sterile Petroleum Gauze, Covidien Petrolatum Gauze, DermaSciences Petrolatum Gauze Non-Adhering Dressing

Reorder No.	Petroleum Gauze Dressing	Quantity	HCPCS
3039	3"×9" – Dressing	4/50/Cs	A6223
3040	3"×9" – Dressing	10/12/Cs	A6223

CuraFoil™ AG Silver Foam Dressing



- Infused with silver ions for antimicrobial protection
- Absorbs up to 5 times its own weight
- Used on partial and full thickness wounds such as pressure ulcers, 1st and 2nd degree burns, surgical wound sites, diabetic ulcers and donor sites

COMPARE TO: DermaRite HydraFoam AG, Smith & Nephew Alleyn AG, Medline Optifoam AG

Reorder No.	CuraFoil™ AG	Quantity	HCPCS
3047	4"×4" – Silver Foam Dressing	10/10/Cs	A4649

DynaFoil™ AG Bordered Silver Foam Dressing



- Self-adherent with silver ions and a secondary adhesive border
- Absorbs up to 5 times its own weight
- Vapor permeable waterproof and leakproof border provides extra protection
- Used on partial and full thickness wounds such as pressure ulcers, 1st and 2nd degree burns, surgical wound sites, diabetic ulcers and donor sites

COMPARE TO: DermaRite HydraFoam AG, Smith & Nephew Alleyn AG, Medline Optifoam AG, Coloplast Biatain AG

Reorder No.	DynaFoil™ AG	Quantity	HCPCS
3049	4"×4"	10/10/Cs	A4649

Xeroform Petrolatum Dressing



- Sterile, fine mesh gauze impregnated with a blend of Xeroform and Petroleum
- Maintains moist healing environment

Reorder No.	Quantity	HCPCS
3053	4"×4" – Dressing	4/25/Cs A6222
3054	5"×9" – Dressing	4/50/Cs A6223

SiliGentle™ Silicone Foam Dressing




- Silicone wound contact layer, absorbent pad layer and a vapor permeable and waterproof film
- Silicone layer removes gently without losing its adherence to minimize skin trauma and pain
- Provides superior absorption and locks exudates away

COMPARE TO: DermaRite ComfortFoam, Coloplast Biatain Silicone Foam, Molnlycke Healthcare Mepilex Border, Smith & Nephew Alleyn Gentle

Reorder No.	SiliGentle™	Quantity	HCPCS
3055	2"×2" – Foam	12/10/Cs	A6209
3056	4"×4" – Foam	12/10/Cs	A6209
3057	6"×6" – Foam	12/10/Cs	A6210
3058	2"×2" – Bordered Foam	12/10/Cs	A6212
3059	4"×4" – Bordered Foam	12/10/Cs	A6212
3060	6"×6" – Bordered Foam	12/10/Cs	A6212
3061	Heel, 10"×9" – Bordered Foam	12/5/Cs	A6213
3062	Sacral, 7"×7" – Bordered Foam	12/5/Cs	A6212
3063	Lg Sacral, 9"×9" – Bordered Foam	12/5/Cs	A6213

DynaGel™ Hydrogel Impregnated Gauze Dressing




- Gentle packing of dry to minimally exuding deep wounds
- Maintains a moist healing environment

COMPARE TO: DermaRite DermaGauze, Skintegrity Hydrogel Impregnated Gauze

Reorder No.	Quantity	HCPCS
3064	2"×2" – Dressing	12/15/Cs A6231
3065	4"×4" – Dressing	12/15/Cs A6231

Oil Emulsion Dressing



- Sterile, knitted gauze dressing impregnated with an oil emulsion blend
- Maintains moist healing environment

Reorder No.	Quantity	HCPCS
3067	3"×3" – Dressing	4/50/Cs A6222
3068	3"×8" – Dressing	6/24/Cs A6223
3070	5"×9" – Dressing	6/12/Cs A6223

Dynarex L-Mesitran Honey Wound Gel



- Creates a moist wound environment
- Facilitates good wound bed preparation

Reorder No.	Quantity	HCPCS
3078	Soft 1.75 oz. (50 g) Tube	24/1/Cs A6250


SiliGentle™ AG Silver Silicone Dressing



- Antimicrobial soft silicone foam absorbs exudate
- Maintains moist healing environment

Reorder No.	Quantity	HCPCS
3080	2"×2" – Non-Adherent	5/10/Cs A6209
3082	4"×5" – Non-Adherent	5/10/Cs A6210
3083	6"×6" – Non-Adherent	5/10/Cs A6210
3084	2"×2" – Bordered	5/10/Cs A6212
3085	4"×4" – Bordered	5/10/Cs A6212
3086	6"×6" – Bordered	5/10/Cs A6213

DynaSorb™ Super Absorbent Dressing



- Super absorbent polymer core rapidly absorbs wound
- Minimizes the risk of maceration and leakage

Reorder No.	Quantity	HCPS
3087	4"×4" – Non-Adherent	12/10/Cs A6196
3088	6"×6" – Non-Adherent	12/10/Cs A6197
3089	4"×4" – Self-Adherent	12/10/Cs A6196
3090	6"×6" – Self-Adherent	12/10/Cs A6197