



Products / Traditional Wound Care / Gauze Dressings / Super Sponges / Sterile Super Sponges Medium



Sterile Super Sponges Medium

2487

5/TY, 120 TY/CS

Product Features

Multi-layer washed, fluff-dried, 100% cotton

Highly absorbent and durable

Thick, diagonally folded construction adds loft to enhance its overall performance

Attributes

| | |
|-------------------|--------------|
| Brand Name | Dukal |
| Size | Medium |
| Ply | 32-Ply |
| HCPC Code | A6402 |
| Sterile | Yes |
| Weight | 11.46757 lbs |

Packaging

| | |
|--------------------------------|--------------------|
| Packaging | 5/TY, 120 TY/CS |
| Inner Packaging | 5/TY |
| Case Packaging | 120 TY/CS |
| Total Case Quantity | 600 |
| Inner Packaging Barcode | (01)10665973002801 |
| Case Packaging Barcode | (01)20665973002808 |

2 Fleetwood Ct.
Ronkonkoma, NY 11779

Phone: 1 (631) 656-3800
Toll Free: 1 (800) 243-0741
Fax: 1 (631) 656-3810

