



Sterile Surgical Drapes, Fenestrated 18" x 26"

20-001

1/PK, 50 PK/BX, 6 BX/CS

Product Features

Made of tissue-poly-tissue

Polyethylene middle layer helps eliminate risk of bacterial migration and moisture strike through

Has a 3" round fenestration

Available in plain for general use and fenestrated to allow access to operative site while protecting the surrounding area

Excellent for minor surgical procedures

Attributes

Brand Name	Dukal
Size	18" x 26"
Sterile	Yes
Weight	17.61493 lbs

Packaging

Packaging	1/PK, 50 PK/BX, 6 BX/CS
Item Packaging	1/PK
Inner Packaging	50 PK/BX
Case Packaging	6 BX/CS
Total Case Quantity	300
Item Packaging Barcode	(01)00665973006734
Inner Packaging Barcode	(01)10665973006731
Case Packaging Barcode	(01)20665973006738