

# Olympia

Portable Table



**STRONGLITE**<sup>TM</sup>  
SERVING MASSAGE PROFESSIONALS SINCE 1986  
Part of the Earthlite family

800.872.0560 • 760.599.1112 • [info@earthlite.com](mailto:info@earthlite.com) • [earthlite.com](http://earthlite.com)



#### Upholstery Color Options


Marie's Beige



Black



#### Adjustable Facecradle



The STRONGLITE™ Olympia Portable Massage Table is made with high quality beech hardwood, coated steel cabling support system and sturdy hardwood corner blocks that makes this a strong and durable treatment table for physical therapy. The Olympia is generously upholstered with 2.5" of small cell foam, and supple soft vinyl with double stitching on its round corners. Solid hardwood beech legs, corner blocks and double knobs plus a strong cabling system and rigid 6mm end panels make it a table with exceptional strength.

It comes with ½-reiki ½-standard endplates that provide the perfect balance between easy access and table stability. A fully adjustable face cradle with a soft, matching face pillow gently supports your clients in the proper massage position. The portable massage table package travels and stores easily in the rugged nylon carry case.

#### Specifications:

- Size: 28" x 73"
- Weight: 39 lbs.
- 2.5" of responsive cushioning
- Adjustable height: 23.5" to 32.5"
- Working weight: 500 lbs.
- 5-year warranty on frame, material and workmanship, 2 years on foam/vinyl

Complete Package includes hardwood beech massage table with round corners, adjustable face cradle, face pillow, matching half round bolster and a rugged nylon carry case with heavy-duty zipper.



**STRONGLITE™**  
SERVING MASSAGE PROFESSIONALS SINCE 1986  
Part of the Earthlite family

800.872.0560 • 760.599.1112 • info@earthlite.com • earthlite.com