*s*tryker[®]

Stryker Prime™ Series

Mobile. Safe. Efficient.



the Stryker Prime Series stretchers are designed to greatly reduce the physical strain on clinicians while enhancing patient care and comfort. Advanced mobility options such as the Zoom® Motorized Drive System or Big Wheel® make transporting heavy patients easy, while electric functionality allows patients to adjust their own position without calling a caregiver for assistance. Fully equipped and highly configurable, the Stryker Prime Series allows you to create the ideal solution to meet your specific needs.



Zoom Motorized Drive System on the Stryker Prime stretcher virtually eliminates manual pushing.



Lift Assist" Backrest uses a drop-seat design which uses the patient's weight to help raise the head of the bed while cradling the patient to help reduce boosting.



Electric functionality includes the ability to give patients control of their own comfort, while Chaperone° Stretcher Exit System helps reduce flight and fall risk. **Advanced mobility,** with optional Big Wheel or Zoom Motorized Drive System, helps make transport safe and efficient.

NEW electric-powered base and litter option makes positioning the stretcher and patient as simple as pressing a button.

Industry's first stretcher exit system alerts the nurse to a potential patient fall.

Optional integrated scale system provides accurate, repeatable patient weights.

Pioneer® mattress offers superior pressure redistribution and comfort as a solution for preventative skin care.

Power-washable* design with 700 lb. capacity.



Stryker Prime Series



Standard Features

- Thermoformed ABS base with extra capacity storage tray
- 700 lb. weight capacity
- Less than 21-inch low height
- Integrated utility tray
- Litter-mounted oxygen bottle holder
- Integrated transfer board
- Power-washable design (1105 and 1115 only)
- Four-wheel, center-locking brake system with dual end activators
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- Glideaway[™] Siderails
- 4 integrated IV receptacles
- +/- Trendelenburg
- 8-inch Omni Surface[™] Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Constant descent hydraulics
- Dual end brake/steer controls
- 2-year warranty

1105: All Standard Features Plus

· Fifth Wheel steering

1115: All Standard Features Plus

• Big Wheel steering

1125: All Standard Features Plus

Zoom Motorized Drive System

1115

Specifications



3800 E. Centre Ave. Portage, MI 49002 U.S.A. t: 269 329 2100 f: 269 329 2311 toll free: 800 787 9537

www.stryker.com

•	26-inch Litter		30-inch Litter

Model Number	1105 (Prime with Fifth Wheel) (€ 1115 (Prime with Big Wheel) (€ 1125 (Prime with Zoom Drive) (€	1105 (Prime with Fifth Wheel)
Overall Length		
(1105 and 1115)	85" (216 cm)	85" (216 cm)
(1125)	86" (218 cm)	86" (218 cm)
Overall Width		
Siderails Up	34" (86.4 cm)	38" (96.5 cm)
Siderails Down (1105 and 1115)	30.25" (76.8 cm)	30.5" (77.5 cm)
Siderails Down (1125)	30.25" (76.8 cm)	31" (78.7 cm)
Patient Surface	26" x 75.25" (66 x 191 cm)	30" x 75.25" (76 x 191 cm)
Weight Capacity	700 lb. (318 kg)	700 lb. (318 kg)
Height Range		
High	34" (86.4 cm)	34" (86.4 cm)
LOW (1105 and 1115)	20.75" (53 cm)	20.75" (53 cm)
Low (1125)	23" (58 cm)	23" (58 cm)
Litter Positioning		
Backrest	0°-90°	0°-90°
Knee Gatch	0°-40°	0°-40°
Trend./Reverse Trend. (1105 and 1115)	±17°	±17°
Trend./Reverse Trend. (1125)	±16°	±16°
Siderails	14" x 57.75" (36 x 147 cm)	14" x 57.75" (36 x 147 cm)
Caster Diameter	8" (20 cm)	8" (20 cm)
EMC	IEC 60601-1-2	

Warranty

Two years parts, labor and travel

Optional Features

- 4-inch (10 cm) Enhanced Comfort Mattress
- 4-inch (10 cm) or 5-inch (13 cm) Ultra Comfort Mattress
- Pioneer Pressure Redistribution Mattress
- Four-sided brake and steer controls (1105 and 1115 only)
- Pop-up push handles (1105 and 1115 only)
- Knee Assist[™] No-crank Hydraulic Knee Gatch
- Exclusive Lift Assist Backrest design with drop seat
- Integrated scale system
- Three-sided hydraulic controls (1105 and 1115 only)

- Defibrillator tray
- Defibrillator tray/foot extender/ charting surface
- Dual-ended siderail release
- Footboard/chartholder
- IV caddy
- IV poles (permanent, folding two- or three-stage)
- Serving tray
- Serving tray holder/footboard
- Siderail pads
- Upright oxygen bottle holder
- Restraint straps
- Retractable cord reel (1125 only)

Electric Functionality

- Electric option (only 30 inches wide)
- Electric patient controls
- Foot end nursing controls
- Powered lift (1105 and 1115 only)
- Chaperone Stretcher Exit System
- 120V or 240V available

Stryker reserves the right to change specifications without notice.

FSC LOGO

Mkt Lit-369 111709 Rev C 2FC/RRP 3M 0112, 11-SPH-034 Copyright © 2012 Stryker