

SUCTION UNITS

Vacuum Levels up to 560 mmHg

Dynarex Resp-O₂™ Suction Units are premium-grade medical devices designed to safely and effectively remove excess secretions and fluids from the airways of respiratory patients. These powerful suction units are easily portable and feature high-quality vacuum regulators and anti-vibration vacuum gauge settings to ensure quiet and efficient operation.

- Designed to provide vacuum levels up to 560 mmHg
- Fiber filter guards against bacteria and contaminants
- Clearly marked gauge allows for easy pressure monitoring
- Oil-less lubrication provides consistent functionality
- Includes an 800 cc disposable suction canister, tubing and bacterial filter



32300



32301



Suction Unit with
Suction Cup Legs

REORDER NO. 32300

PACKAGE QUANTITY
1/Case

Suction Unit
with Base

REORDER NO. 32301

PACKAGE QUANTITY
1/Case

