

BRODA[®]

Synthesis

Positioning Wheelchair



Shown with rear 14" mag wheel option

Patient, Facility, & Caregiver Benefits

- Increases quality of life and engagement with surroundings
- Promotes safe patient handling with enhanced positioning
- Proper seating alignment reduces patient fatigue and discomfort



Tilt-in-space functionality provides additional pressure offloading



Front pivot seat tilt maintains a forward line of sight, enhancing engagement.



Adjustable upper lateral supports fit a variety of users



Optional WC-19 Transportation Package enables safe patient travel

brodaseating.com | 1-844-552-7632

A FULL RANGE OF POSITIONING FEATURES
WITH CONTEMPORARY DESIGN

UPHOLSTERY OPTIONS



SPECIFICATIONS

| | |
|-----------------------------------|--|
| Tilt - Infinitely Adjustable | up to 40° |
| Recline - Infinitely Adjustable | up to 90° |
| Weight Capacity | 350 lbs. (158 kg) |
| Seat Widths | 16", 18", 20", 22", 24", 26" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm) |
| Seat Depths | 17", 18.5", 20", 22" (43cm, 46cm, 50cm, 56cm) |
| Floor-to-Seat Heights | 19", 20.5", 22", 23.5" (48cm, 51cm, 55cm, 59cm) |
| Wheels | 14", 20" Rear Mount Mag (35cm, 50cm, 55cm) 14" Mid Mount Mag (35cm) |
| Lifetime Warranty Frame | 10 years on 16 Gauge tubular steel frame |
| Pressure-Reducing Seat & Backrest | Integrated with Comfort Tension Seating® |
| Standard Padding | Neckrest, upper laterals, side pads, sole and calf pads, seat and back pads |

FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning straps