

Large Top Bariatric Step Stool



MODEL: T-6242


- Durable chrome finish
- Extra wide and deep top
- Welded reinforced "X" brace
- 600 lbs. load capacity under normal use

Shipping weight: 14 lbs.

SPECIFICATIONS

PLATFORM: 17.25" L X 12.25" W X 9"H

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1