


PARALLEL BAR, BARIATRIC



IEC60601-1(ed.3)
RM 14971



1000 lb capacity

Tapered Extension

<u>Item No.</u>	<u>Length</u>	<u>Platform Width</u>
TWG0542.BA	1 ft	54"W

FEATURES:

- 1000 lb Weight Capacity
- 1 ft Tapered Platform Extension
- Steel Frame

Base: 10 Gauge Steel; 1 3/4" high
Finish: Pearl Grey Standard; Porcelain Optional.
Warranty: Ten (10) Year Limited Warranty.

*Nominal Dimensions Stated