



Deluxe Sphygmomanometer, Black

2010

1/EA

Product Features



WARNING: This product can expose you to chemicals including lead, which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.

Measures blood pressure from 0mm to 300mm Hg

Black nylon cuff with index and range markings

Two tube non-latex bladder with deluxe inflation system features a large knurled air release valve

Both the air release valve and end valve are filtered to keep the unit free from dust and dirt

Black nylon carry case

Cuff dimensions: 20.5" W x 5.5" H

Attributes

Brand Name	Tech-Med
Size	Adult
Color	Black
Material	Nylon
Sterile	No
Weight	1.05639 lbs

Packaging

Packaging	1/EA
Item Packaging	1/EA
Inner Packaging	12EA/BX
Case Packaging	2BX/CS
Total Case Quantity	24
Item Packaging	(01)00665973016917