

CARPAL TUNNEL GLOVE



THERMOSKIN

TRIOXON[®] ADVANTAGE

Exclusive Trioxon 3 dimensional lining provides heat therapy & pain relief. Soft, comfortable inner layer wicks away moisture. Middle layer helps increase skin temperature for heat therapy. Outer layer provides light compression.



Carpal Tunnel Glove combines heat therapy and light compression for arthritis or aching hands. A metal splint positioned at 25 degrees provides the protection and support needed for Carpal Tunnel Syndrome, overused or weak wrists. The outer layer of fabric is textured to provide additional grip to aid in daily living activities. A fully adjustable locking strap wraps around the wrist for a better fit and to keep the glove in proper position.

PRODUCT ADVANTAGES:

- Lined with the exclusive Trioxon[®] Advantage material for heat therapy & enhanced recovery
- Clinically tested to increase skin temperature and blood flow for pain relief
- Provides light but firm compression to counteract tissue swelling
- Metal splint at the palmer (inner) surface of the wrist for stability
- Combines heat therapy for pain relief with a Carpal Tunnel brace.

SIZING GUIDE

Measure circumference around knuckles



PART NUMBERS		SIZE	INCHES
LEFT	RIGHT		
82197	82198	XS	6-6¼"
83197	83198	S	7-7¼"
84197	84198	M	8-9"
85197	85198	L	9¼-10½"
86197	86198	XL	10¾-11½"
87197	87198	XXL	11¾"+

THERMOSKIN TECHNOLOGY

Thermoskin provides effective treatment and pain relief for Sports Injuries, Arthritis, Repetitive Strain Injuries and many other conditions. Thermoskin combines compression, support and naturally generated body heat for therapy and pain relief.

The benefits of Thermoskin Technology include:

- Promotion of increased blood flow facilitates enhanced recovery.
- Provides light but firm compression to counteract tissue swelling.
- Clinically tested to increase the skin and muscle temperature.
- Provides heat therapy and pain relief.
- Increases elasticity and reactivity of the muscle.

HEAT THERAPY

Thermoskin is clinically proven to increase skin and muscle temperature. Providing heat to the affected area is one of the best modalities employed in the treatment of soft tissue injuries. Heat also provides effective relief for the pain and stiffness associated with arthritis.

TRIOXON® ADVANTAGE

Exclusive 3 dimensional lining that provides insulation and wicking of moisture via air circulation. This allows the skin to remain well oxygenated and comfortable for extended periods.



CLINICALLY PROVEN


An independent University study on the effect of Thermoskin devices on local blood flow, muscle function and microclimate produced these key findings:

- Wearing a Thermoskin Knee product **increased blood flow** by 8% and wearing a Thermoskin Ankle/Foot product increased blood flow by 19%.

Note: Increased blood flow moves oxygen rich blood to the area and facilitates the clearing of cellular debris, thereby helping with the healing process.

- Wearing various Thermoskin products such as a Knee or Lumbar Support provided substantially **greater skin temperature** when compared to the control condition.

Note: Elevated skin temperature works as heat therapy thereby assisting in the healing process.



Thermal image of right foot showing skin temperature.

- Wearing a Thermoskin Knee product provided **greater muscle temperature** when compared to the control condition.
- Wearing Thermoskin Shorts **improved flexibility** at the hip flexors when compared to the control condition.

Note: Increased muscle temperature promotes optimal muscle function by increasing the elasticity and reactivity of the muscle.

- Wearing a Thermoskin Wrist product generated less relative humidity than our direct competitor.

Note: The greater the Relative Humidity the more we sweat. Thermoskin wicks away moisture and allows for a free flow of air between the material and the skin. This allows the skin to remain well oxygenated and comfortable for extended periods.