

Three Piece Transfer Tub Bench

HCPCS: E0247-1



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Tool-free assembly of back, legs, and arm (Figure A - C)
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

PRODUCTS	Description
RTL12031KDR	Tranf Bnch, 3 pc.Assembly,RTL,1 EA
SPECIFICATIONS	Description
Outside Legs	32" (W) x 23" (D)
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Seat Height	17.5" - 22.5"
Shipping Weight	14.37 lbs
Warranty	Limited Lifetime
Weight	12.48 lbs
Weight Capacity	400 lbs