

Guide to 3pp® Thumb Splints Web Pages

Everything you need to know to about 3-Point Products Thumb Splints can be found in the following web pages on our website (www.3pointproducts.com), including product details, testimonials, instructional videos, FAQ's, blog posts and where to buy information.

3pp® ThumSling®

www.3pointproducts.com/3pp-thum-sling/



Wrap around strap support the CMC joint to reduce subluxation and pain.

3pp® ThumSling® Long

www.3pointproducts.com/3pp-thumsling-long/



Wrist and CMC support provides comfort and function.

3pp® ThumSling® NP

www.3pointproducts.com/3pp-thumsling--np/



Stretch resistant material and a longer design provides support to reduce CMC joint pain.

3pp® ThumSling® Long NP

www.3pointproducts.com/3pp-thumsling--long-np-0/



The next level of control for wrist and CMC support.

3pp® ThumSpica™

www.3pointproducts.com/thum-spica/



One-size wrap-on support positions the thumb for rest or function.

3pp® ThumSpica™ Plus

www.3pointproducts.com/3pp-thumspica-plus/



Combines a custom molded insert and easy-to-fit wrap-on support.

ThumSaver™ MP

www.3pointproducts.com/thum-saver-mp/



Stabilizes the MP and controls CMC motion.

ThumSaver™ CMC Short

www.3pointproducts.com/cmc-short/



Provided firm thumb support and limits wrist motion.

ThumSaver™ CMC Long

www.3pointproducts.com/cmc-long/



Carefully trimmed to leave the ulnar side of the hand free.

3pp® Ultra Spica

www.3pointproducts.com/ultra-spica/



Controls MP motion and supports the CMC.

3pp® ThumDuction™ Strap

www.3pointproducts.com/thum-duction-strap/



Two-part straps that position the thumb for function and control of muscle tone.

3pp[®] Thumb Splints Sizing & Ordering Information

3pp[®] ThumSling[®]

www.3pointproducts.com/3pp-thum-sling/

Lay hand flat and measure around the MP joints (large knuckles).

Size	Color	Right	Left	Fits (Circumference at MP's)
Sm/Med	Black	P2006-R23BK	P2006-L23BK	6" to 7 1/4" (15.25-18.4cm)
Sm/Med	Grey	P2006-R23GR	P2006-L23GR	6" to 7 1/4" (15.25-18.4cm)
Med/Lrg	Black	P2006-R34BK	P2006-L34BK	7 1/4" to 9 1/4" (18.4-22.8cm)
Med/Lrg	Grey	P2006-R34GR	P2006-L34GR	7 1/4" to 9 1/4" (18.4-22.8cm)

3pp[®] ThumSling[®] NP

www.3pointproducts.com/3pp-thumsling--np/

Lay hand flat and measure around the MP joints (large knuckles).

Size	Right	Left	Fits
Sm/Med	P2010-R23	P2010-L23	6" to 7 1/4" (15.25-18.4cm)
Med/Lrg	P2010-R34	P2010-L34	7 1/4" to 9 1/4" (18.4-22.8cm)

3pp[®] ThumSpica[™]

www.3pointproducts.com/thum-spica/

Lay hand flat and measure around the MP joints (large knuckles).

Size	Item Number
One Size	P3009-R (Right) Fits 6"-9" MP Circumference (15.25-22.8cm)
One Size	P3009-L (Left) Fits 6"-9" MP Circumference (15.25-22.8cm)

ThumSaver[™] MP

www.3pointproducts.com/thum-saver-mp/

Measure the length from the IP joint (knuckle at the end of the thumbs) to the wrist. Measure around the IP joint.

Size	Right	Left	Length	Circumference
Sm	P2001-R2	P2001-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P2001-R3	P2001-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P2001-R4	P2001-L4	4" to 4 1/2"	3 1/8" to 3 1/4"
X-Lrg	P2001-R5	P2001-L5	4 1/2" to 5"	3 1/4" to 3 1/2"

ThumSaver[™] CMC Long

www.3pointproducts.com/cmc-long/

Measure the length from the IP joint (knuckle at the end of the thumbs) to the wrist. Measure around the IP joint.

Size	Right	Left	Length	Circumference
Sm	P3000-R2	P3000-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P3000-R3	P3000-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P3000-R4	P3000-L4	4" to 4 1/2"	3 1/8" to 3 1/4"

3pp[®] ThumDuction[™] Strap

www.3pointproducts.com/thum-duction-strap/

Measure length from the end of the joint of the thumb to the wrist. Fits both right and left hands.

Size	Item Number	Fits
Adult	P2005	2 3/4" to 4 1/2"

3pp[®] ThumSling[®] Long

www.3pointproducts.com/3pp-thumsling-long/

Lay hand flat and measure around the MP joints (large knuckles).

Size	Item Number
One Size	P3026-R (Right)
One Size	P3026-L (Left)

3pp[®] ThumSling[®] Long NP

www.3pointproducts.com/3pp-thumsling--long-np-0/

Lay hand flat and measure around the MP joints (large knuckles). One size fits most.

Size	Right	Left	Fits
One Size	P3010-R	P3010-L	6" to 9" MP Circumference(15.25-22.8cm)

3pp[®] ThumSpica[™] Plus

www.3pointproducts.com/3pp-thumspica-plus/

One size fits most. Measure around large knuckles of hand.

Size	Item Number	Fits
One Size	P3019-R (Right)	6" to 9" MP Circumference(15.25-22.8cm)
One Size	P3019-L (Left)	6" to 9" MP Circumference(15.25-22.8cm)

ThumSaver[™] CMC Short

www.3pointproducts.com/cmc-short/

Measure the length from the IP joint (knuckle at the end of the thumbs) to the wrist. Measure around the IP joint.

Size	Right	Left	Length	Circumference
Sm	P2000-R2	P2000-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P2000-R3	P2000-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P2000-R4	P2000-L4	4" to 4 1/2"	3 1/8" to 3 1/4"
X-Lrg	P2000-R5	P2000-L5	4 1/2" to 5"	3 1/4" to 3 1/2"

3pp[®] Ultra Spica

www.3pointproducts.com/ultra-spica/

Measure the length from the IP joint (knuckle at the end of the thumbs) to the wrist. Measure around the IP joint.

Size	Right	Left	Fits
Sm	P3100-R2	P3100-L2	3" to 3 1/2"
Med	P3100-R3	P3100-L3	3 1/2" to 4"
Lrg	P3100-R4	P3100-L4	4" to 4 1/2"