

Prevent Falls & Wandering

Timed Pads for Bed & Chair Exit Alarm



Works with:

Change Pad Indicator
Fall Monitor
TL-2100CP

Eliminate F-Tags!

Eliminate F-Tags when using timed pressure pads. When used with the Change Pad Indicator Monitor, automatically tracks the life of the sensor pad and alerts caregiver 10 days prior to expiration date.

Timed Chair Pads

- Unique coiled cord reduces tripping and entanglement hazards

Timed Bed Pads

- Unique Safe-T-Release cord separates instead of breaking or pulling away from the pad

MFR#	SMART PRODUCT DESCRIPTION
TB-WI	Timed Bed Pad 20"x30" (1 year warranty)
TB-90W	Timed Bed Pad 20"x30" (90 day warranty)
TB-45W	Timed Bed Pad 20"x30" (45 day warranty)
TB-RI	Timed Bed Pad 10"x30" (1 year warranty)
TB-90R	Timed Bed Pad 10"x30" (90 day warranty)
TB-45R	Timed Bed Pad 10"x30" (45 day warranty)
TC-WI	Timed Chair Pad 10"x15" (1 year warranty)
TC-90	Timed Chair Pad 10"x15" (90 day warranty)
TC-45	Timed Chair Pad 10"x15" (45 day warranty)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com