

TranSport Aluminum Transport Chair

HCPCS: E1038



FIGURE A



FIGURE B

- Weighs 17.5 lbs without footrests
- Aluminum frame is lightweight and strong
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried
- Padded flip-back arms are removable
- Upholstery is comfortable, lightweight and easy-to-clean
- Two-way back release allows to push or pull lever
- The transport chair conveniently folds for travel and storage (Figure A)
- Back folds forward
- Comes with tool-free height adjustable swing-away footrests
- Aluminum rear wheel locks
- Weight Capacity: 250 lbs

PRODUCTS

Description

TS19	Tran Sport Chair,FlipBack-RemArms,1/cs
------	--

SPECIFICATIONS

Description

Back Height	16"
Carton Shipping Weight	19.5 lbs
Overall Height	36"
Overall Length w/ Riggings	38"
Overall Width	24.5" (W) (open) 10" (W) (closed)
Seat Dimensions	19" (W) x 16" (D)
Seat To Floor Height	18"
Warranty	Limited Lifetime
Weight Capacity	250 lbs
Weight Without Footrests With Footrests	17.5 lbs 20 lbs

PARTS

Description

STDS871	Hand Grip,Sentra,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
STDS2A1800	Rear Wheel, DFL, TS19,1/ea, EA
STDS2A1801	REAR HUB CAP,DFL,1 PR, PR
STDS2A184L	BRAKE ASSY LEFT,DFL,1/EA, EA
STDS2A1811	REAR AXLE BOLT,DFL,1 PR, PR
STDS2A1806	FRONT AXLE BOLT,DFL,1 PR, PR

PARTS	Description
STDS2A1810	FRONT WHEEL,DFL,TS19,1/ea, EA
TRB-1	BackSupportReleaseLever./Trans,1/ea, EA
STDS1A2307	Seat Upholstery,TS19,1/ea, EA
STDS2A184R	BRAKE ASSY RIGHT,DFL,1/EA, EA
STDS1A2308	Back Upholstery,TS19,1/ea, EA
STDS001013	Folding Cane Hindge TR18,TS19,1/ea, EA
STDS1A2305R	RightBackCaneAssy,TS19,1/ea, EA
STDS2A1803	Cross Frame Assist DFL,1/EA, EA
STDS2A1807	SEAT RAIL PAD,DFL,NITRO,1/EA, EA
STDS1A2318R	Brake Assembly, Right, TS19, 1/ea, EA
STDS1A2309R	RightFootrestAssembly,TS19,1/ea, EA
STDS1069-4	Hanger Pins, 4 per PKG, PKG