




TUBES AND CONNECTORS

Dynarex Resp-O₂™ Flex Tube/Flex Connectors add length to vertical circuits and HMEs to help improve system flexibility and range of motion.

- Provides added flexibility to ventilator circuits
- Lightweight, pliable design
- Not made with natural rubber latex or DEHP.



Universal 6"
Flex Tube

REORDER NO. 30600

DEAD SPACE
40 mL

WEIGHT
9.5 g

CONNECTIONS
15F/22F

PACKAGE QUANTITY
50/Case



Tracheostomy
Adult Flex Connector

**15 MM OD × 22 MM
OD/15 MM ID**

REORDER NO. 30601

DEAD SPACE
8 mL

WEIGHT
6.5 g

CONNECTIONS
22M/15F-15M

PACKAGE QUANTITY
50/Case



Tracheostomy
Pediatric Flex
Connector

15 MM OD × 15 MM ID

REORDER NO. 30602

DEAD SPACE
8 mL

WEIGHT
5.5 g

CONNECTIONS
15F-15M

PACKAGE QUANTITY
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.



©2021 Dynarex Corporation

R210604