



Universal 5" Walker Wheels

HCPCS: A9900 OR E0155



- Converts folding walker into wheeled walker
- Fits most manufacturers walkers
- Allows for 8 height adjustments
- Rubber wheel allows walker to roll easily and smoothly over irregular surfaces
 Comes with rear glide caps, (item #10107) and glide covers (item #10107C) allowing use on all surfaces
- 350 lb Weight Limit

PRODUCTS	Description
10109	Walker Wheels,Univrsl 5",1pr/bx, RTL
SPECIFICATIONS	Description
Weight Capacity	350 lbs