

XL Padded Blood-Drawing Chair

2574/2575



ISO 9001-2008 Certified

A strong, sturdy bariatric chair with all the features and quality you expect.

Steel constructed frame provides strength and durability for bariatric use. The padded seat and back give exceptional patient comfort. Winco's commitment to quality will provide you with years of trouble-free service.



2575



2574

KEY FEATURES

- Padded seat and back
- Contoured urethane arms
- Sturdy leveling feet
- Adjustable height armrests
- Laminated cabinet (2574 only)
- Drawer roller glides (2574 only)

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

XL PADDED BLOOD-DRAWING CHAIR 2574/2575



- 1 Comfortably padded seat and backrest.
- 2 Flip-up, contoured urethane armrest provides easy access and added safety for the user.
- 3 Adjustable height armrests.
- 4 High impact leveling feet ensure a sturdy chair on almost any surface.
- 5 Powder-coated steel frame.


SPECIFICATIONS

2574/2575

Seat Height	22.5" (57.2cm)
Seat Width	24.25" (61.6cm)
Seat Depth	18.5" (47cm)
Overall Height	36.5" (92.7cm)
Overall Width	2574: 41.5" (105.4cm)
	2575: 34" (86.4cm)
Back Height Above Seat16" (40.6cm)
Overall Length25" (63.5cm)
Adjustable Arm Height Range	28.25" - 34.25" (71.76 - 87cm)
Weight	2574: 72 lbs. (32.7kg)
	2575: 55 lbs. (24.9kg)
Weight Capacity	450 lbs. (204kg)

*Weight shown as a minimum-maximum range & varies depending on chair options ordered.

OPTIONS



TB-133 **TB**



Pivot Arm **PR**
(both models)
Dual Pivot Arm **P2**
(2575 only)



IV Pole **IV**
(right-rear only)

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS

01-26-11

REV.B

009055