



*The Winco Privacy Screen Collection*

Privess™



Privess™ is a broad range of patient privacy screens providing flexible solutions to meet the wide spectrum of requirements for healthcare professionals, military procurement, international aid or disaster relief coordinators and educational institutions worldwide.

- Cost-effective selections, from innovative modular to attractive designer models
- Permanently mounted, stationary, portable and transportable configurations available
- Sturdy steel or lightweight aluminum construction options with easy-clean smooth panels

Privess™ represents our ongoing dedication to privacy so clinicians can confidently focus on their patients. For every conceivable need, budget or profession, Privess™ has a custom tailored answer to your patient privacy needs.




## **PRIVESS™ BASIC**

### **3-PANEL PRIVACY SCREEN | Model 3130**

For healthcare, disaster relief and education professionals who require a sturdy metal screen that can easily be moved to protect patient privacy. The Privess™ Basic is a cost effective durable steel 3-panel screen that features:

- Casters for mobility and crutch tips for stationary use
- Smooth white vinyl or optional Sure-Chek® panels for easier cleaning
- Durable steel 3-panel design

Unlike similar style aluminum competitors the Privess™ Basic 3-Panel Privacy screen offers portable, high-quality steel framed privacy protection in a value priced package



## **PRIVESS™ MODULAR**

### **PRIVACY SCREEN | Model 3730 / 3740**

Value-minded healthcare professionals, military, disaster relief coordinators and educational institutions will appreciate these cost effective, flexible, transportable, easy-to-clean solutions for patient privacy. The Privess™ Modular screens feature:

- Unlimited configurations within its 3 and 4 panel models.
- Easy to clean smooth vinyl panels
- Standard white vinyl or optional Sure-Chek® colors available

The Privess™ Modular privacy screen's unique snap together design offers 360 degree swivel and infinite connectivity. Its compact storage footprint, lightweight and flexible design is ideal for rapid deployment, temporary use situations.



## **PRIVESS™ UNIPANEL**

### **PRIVACY SCREEN | Model 3420**

Perfect for healthcare professional, disaster relief and educational professionals who work within limited confines and desire a compact portable solution to protect patient privacy. The Privess™ Unipanel is a compact robust single panel screen that features:

- Casters for mobility or crutch tips for stationary use
- Smooth white vinyl or optional Sure-Chek® panels for easier cleaning
- Robust, durable steel construction

The Privess™ Unipanel is a right-sized single panel solution perfect for patient privacy in smaller areas or facilities with little storage room. Privess™ Unipanel is a quality steel framed alternative to inexpensive, lightweight, imported aluminum solutions.



## **PRIVESS™ SWING-AWAY**

### **PRIVACY SCREEN | Model 3400**

For healthcare and educational professionals seeking a fixed placement privacy solution that can easily be drawn to protect patient privacy or quickly swung out of the way for unobstructed viewing and patient access. The Privess™ Swing-Away privacy screen is a wall mounted, telescoping privacy screen that:

- Provides 180 degree range of motion (swivel)
- Telescopes from 40" to 92.5"
- Standard white vinyl or optional Sure-Chek® colors available

The Privess™ Swing-Away privacy screen is convenient to use and stores completely out of the way yet provides a wall of privacy on demand. Unlike ceiling mounted track based screens, this Swing-away model will not interfere with patient lifts or other equipment.



## **PRIVESS™ ELITE**

### **PREMIUM PRIVACY SCREEN | Model 3170**

For quality and style minded users who need attractive, décor matched, stable fixed or portable quality patient privacy solutions. Privess™ Elite is a high quality screen featuring:

- Durable steel construction
- Quality designer grade vinyl colors to enhance any décor
- Casters for mobility or crutch tips for stationary use

Privess™ Elite construction showcases significantly upgraded durability and flexibility unlike imported unstable lightweight construction with only a thin sheet vinyl. The Privess™ Elite offers a quality designer product that is at home in healthcare institutions and surgery centers, physician offices, resort style spas and even retail settings where stylish substantial privacy screens are desired.



# The Winco Privacy Screen Collection

## Privess™


### ANATOMY OF THE PRIVESS™ PRIVACY SCREEN

#### Includes Models 3130 / 3170\* / 3420


\* 3170 available in all Winco standard vinyl colors or Sure-Chek®



#### Includes Model 3400



#### Includes Models 3730 / 3740



#### SURE-CHEK®

All Privess™ screens are available with optional Sure-Chek® anti-microbial vinyl. Beyond the important anti-microbial features, Sure Chek® is also anti-static, odor, stain and flame resistant and is available in 5 decorative colors.

## MODELS & SPECIFICATIONS

Available in a variety of colors



MODEL	3130	3420	3400	3730 / 3740	3170
Base Size	20" x 27"	20" x 42"	N/A	15" x 40" / 15" x 70"	23" x 31"
Overall Height	70"	69"	72" Recommended	67.5" / 67.5"	70"
Overall Width	79"	42"	40" – 92.5"	60.5" / 80.5"	84"
Single Panel Length	49"	49.5"	69"	45.5" / 45.5"	55"
Single Panel Width	27"	42"	92.5"	18.5" Vinyl, 20" Frame	28"
Weight	25 lbs	15 lbs	8 lbs	7 lbs / 9.25 lbs	45 lbs

## FEATURES

MODEL	PORTABLE	WALL MOUNT	MODULAR	SWING / TELESCOPING	STEEL	ALUMINUM
3130	Y	N/A	N	N	Y	N
3170	Y	N/A	N	N	Y	N
3420	Y	N/A	N	N	Y	N
3400	N	Y	N	Y	Y*	Y*
3730	Y	N/A	Y	N	N	Y
3740	Y	N/A	Y	N	N	Y

\* Wall bracket and support are steel, extension rod is aluminum.

## OPTIONS

CRUTCH TIP	CASTERS	WHITE COLOR	SURE-CHEK® COLORS	VINYL COLORS
S	S	S	0	N/A
S	S	S	0	S/0
S	S	S	0	N/A
N/A	N/A	S	0	N/A
S	N/A	S	0	N/A
S	N/A	S	0	N/A

S-Standard 0-Optional N/A-Not Applicable

S/0-Grade A vinyl colors are standard. Grades B, C and Silvertex are optional

## RECOMMENDED ENVIRONMENTS

MODEL	WHERE PORTABILITY IS CRITICAL	FIXED POSITION, YET MOVES OUT OF THE WAY	WHEN WHEELS ARE NEEDED	WHERE LIGHTWEIGHT IS IMPORTANT	WHERE MULTIPLE SCREENS ARE NEEDED	FOR CUSTOM SHAPES (CIRCLE, ROOM, ETC)	RAPIDLY DEPLOYED
3130	○		●				
3170	○		●				
3420	○		●				
3400		○					○
3730	●			●	●	●	
3740	●			●	●	●	

● Optimal Application ○ Suitable Application

\*Winco Privacy Screens are not intended for wet environments and should never be used as shower curtains. Winco Privacy Screens are not recommended for outdoor use.

Sure-Chek® is a registered trademark of Herculite Products Inc.



## Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com