

ECG Electrodes

Ag/AgCl Pregelled Monitoring Electrodes

Contains
4

AW-R2246-04 Rev. A



ZOLL



ZOLL Medical Corporation
269 Mill Road
Chelmsford, MA 01824-4105 USA
800-348-9011 Made in USA



ZOLL International Holding B.V.
Newtonweg 1B 6662 PV ELST
The Netherlands Tel. +31 (0) 481 366410



2

Cautions:

- Do not use if gel is dry.
- Replacement of electrodes after each 24 hour period is suggested. Vary placement to minimize skin sensitivity.



3



4

