

**Zoo Bus, Pediatric Ready Room**



*an all-inclusive Clinton exclusive*

Who doesn't like friendly animals? And better yet, who doesn't like the convenience of an entire pediatric exam room designed for you? Get some of the best young patient and parent pleasing Clinton Kids<sup>™</sup> products all packaged at a price that is ready to go.

*(no model substitutions)*

## MODEL 7020-RR

---

### PACKAGE INCLUDES:

- 7020 Fun Series Zoo Bus Treatment Table with Black upholstery
- 6120-BW Schoolhouse Base & Wall Cabinets with Slate Gray top and sides plus stainless steel sink and chrome faucet set (choose right or left location or no sink. If no location is selected, sink will be placed on left as shown.)
- 10-K Koala Pediatric Laminate Step Stool
- GS-3022 Double, Slimline, Stainless Steel,\* Glove Box Holder
- 21335 in your choice of upholstery colors
- C-40B Premium, Padded, Side Chair in your choice of upholstery colors

\* Not MRI compatible

### 7020 Zoo Bus Treatment Table

---

- Clinton original illustrations
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on two long side and two short ends, manufactured as pictured
- Optional locking doors available
- All colors are as shown
- 3 Storage compartments with one adjustable shelf each
- Soft-close, door hinges
- Chip resistant PVC edges
- Seamless rounded corner top
- Premium, stain resistant, woven knit back upholstery (black)
- Adjustable feet
- Doors are only on one side
- Paper dispenser included
- Special order reverse side doors are available
- 350 lbs. load capacity under normal use

#### SPECIFICATIONS

72"L x 27"W x 34"H



#### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov). **1**

### 10-K Koala, Yellow, Laminate Step Stool

---

- Heavy laminated wood construction
- Non-skid tread
- Coordinates with Zoo Bus & Kangaroo Country tables
- Koala graphic on both long sides

#### SPECIFICATIONS

20"L x 16"D x 7"H



#### Prop 65 Warning

**!** **WARNING:** This product can expose you to chemicals, including vinyl chloride, which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov). **3**

## 6120-BW Schoolhouse Cabinets

### "Schoolhouse" Base Cabinets

4 doors with adjustable shelves and 3 drawers.  
Shipped in one unit.

### "Schoolhouse" Wall Cabinets

4 doors with adjustable shelves. Shipped in one unit.

### Base & wall cabinets are only sold in sets

- Clinton original illustrations
- Pre-assembled units save time
- Easy-clean all laminate construction
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Tough, chip-resistant PVC edges
- Assembled with steel fasteners
- 4 adjustable leveling feet on each unit
- 3-way adjustable, soft-close 110° hinges
- Flat panel drawer fronts with metal sides that hold up to 75 lbs.
- Recessed toe kick
- Easy hook and rail installation for upper cabinets
- Slate Gray edgebanding on all edges
- Black chrome door/drawer pulls
- Customizable with 4 options
- Easily removes for relocation
- Drawers can be added to any base unit
- Doors and drawer can be moved to accommodate sink location
- Shown with optional sink and faucet

#### SPECIFICATIONS

BASE: 61"L x 18"D x 35"H

WALL: 61"L x 12"D x 24"H



### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

## GS-3020 Double, Stainless Steel,\* Glove Box Holder

- Stainless steel construction
- Open design allows for easy glove identification
- All units mount horizontal or vertical
- Mounting hardware included
- Universal size fits most glove boxes

*\*Clinton stainless steel items are not M.R.I. compatible.*



Shown in the horizontal position

#### SPECIFICATIONS/Horizontal

10.5" L x 3.75" D x 10" H

### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



## 2135 Standard Series, 5-Leg, Pneumatic Stool

- Thick padded seat for comfort
- 24" Black nylon, 5-leg base
- Nylon dual wheel casters
- Easy, pneumatic height adjustment
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- Maximum load capacity 250 lbs. under normal use



### Prop 65 Warning

**WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

2

#### SPECIFICATIONS

Seat Diameter	16"
Base Diameter	24"
Height Range	18"-23"

## C-40B Premium, Padded, Black Frame Side Chair Without Arms

- Black powder coated frame finish
- Heavy duty 1-1/2" diameter legs
- Solid plywood seat and back
- Wall saver guard
- 2" thick, firm foam padded seat and back
- Contoured seat and back
- Adjustable leg levellers
- Premium, knit-backed upholstery
- 300 lbs. load capacity under normal use



### Prop 65 Warning

**WARNING:** This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: [www.P65Warning.ca.gov](http://www.P65Warning.ca.gov).

7

#### SPECIFICATIONS

Seat	18" W x 18" D
Seat Height	18"

#### Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black